

Product datasheet for **SC206455**

GSTM2 (NM_000848) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: GSTM2 (NM_000848) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: GSTM2
Synonyms: GST4; GSTM; GSTM2-2; GTHMUS
ACCN: NM_000848
Insert Size: 480 bp
Insert Sequence: >SC206455 3'UTR clone of NM_000848
The sequence shown below is from the reference sequence of NM_000848. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG  
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC  
ACAAAGATGGCTGTCTGGGCAACAAGTAGGGCCTTGAAGGCCAGGAGGTGGGAGTGAGGAGCCCATAC  
TCAGCCTGCTGCCAGGCTGTGCAGCGCAGCTGGACTCTGCATCCCAGCACCTGCCTCCTCGTTCCTTT  
CTCCTGTTTATTCCCATCTTTACTCCCAAGACTTCATTGTCCCTTTCACTCCCCCTAAACCCTGTCC  
CATGCAGGCCCTTTGAAGCCTCAGCTACCCACTATCCTTCGTGAACATCCCCTCCCATCATTACCCTTC  
CCTGCACTAAAGCCAGCCTGACCTTCCTTCTGTAGTGGTTGTGTCTGCTTTAAAGGCCCTGCCTGGC  
CCCTCGCCTGTGGAGCTCAGCCCCGAGCTGTCCCCGTGTTGCATGAAGGAGCAGCATTGACTGGTTTAC  
AGGCCCTGCTCCTGCAGCATGGTCCCTGCCTTAGGCCTACCTGATGGAAGTAAAGCCTCAACCACA  
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA  
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_000848.4](#)



[View online »](#)

Summary: Cytosolic and membrane-bound forms of glutathione S-transferase are encoded by two distinct supergene families. At present, eight distinct classes of the soluble cytoplasmic mammalian glutathione S-transferases have been identified: alpha, kappa, mu, omega, pi, sigma, theta and zeta. This gene encodes a glutathione S-transferase that belongs to the mu class. The mu class of enzymes functions in the detoxification of electrophilic compounds, including carcinogens, therapeutic drugs, environmental toxins and products of oxidative stress, by conjugation with glutathione. The genes encoding the mu class of enzymes are organized in a gene cluster on chromosome 1p13.3 and are known to be highly polymorphic. These genetic variations can change an individual's susceptibility to carcinogens and toxins as well as affect the toxicity and efficacy of certain drugs. [provided by RefSeq, Jul 2008]

Locus ID: 2946

MW: 17