

Product datasheet for **SC206545**

C15orf40 (NM_001160114) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: C15orf40 (NM_001160114) Human 3' UTR Clone
Symbol: C15orf40
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_001160114
Insert Size: 514 bp
Insert Sequence: >SC206545 3'UTR clone of NM_001160114
The sequence shown below is from the reference sequence of NM_001160114. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTTGTTGAATGCAGTTTTCTAGATTGTAATTAATAATCATATCTCCACAACGGTTCTGGAATTCTGG
GGATTGCCTTCTAATCCAGTTTACAAGTGACAGAAAACCTGTTATGCTGCTTTGTAGTCTTACCTTAGA
TTTTTAAATGCTTATTTAAATGTAATAAATTTTTTTTTTACATCTTTAAAGACTATATTTACCATT
TTGTTGAAAAAGTAAAGAAAAATCATGACTTGTAATTAATGTGAATTATAATTTTTGAGAAGTCACTG
AATTTTCATCTTCTCTAATAACTTAAATAATCAGCAAATTCAGTAGTTTTTATGAATCTCCAAAGATC
CTTAAAAGAAGAAAGAGGACTTGGAAGCTATTTTCAGTATTTGTAAGCTTGTCAGAAATTACATCTA
TATTAAGTTTCAAAGCTAATAATATATTTTTACATCTTCAGGTCAAGTTATAGTAATTTATTAAGA
AAATAAATTAATTTCTTTGCTGATCAGGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online >](#)

RefSeq: [NM_001160114.2](#)

Locus ID: 123207

MW: 20.4