

Product datasheet for **SC206571**

C2orf42 (NM_017880) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: C2orf42

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_017880

Insert Size: 501 bp

Insert Sequence: >SC206571 3'UTR clone of NM_017880
The sequence shown below is from the reference sequence of NM_017880. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCCTCTGACCACTATTACTTTCCCTAAAGCAAAAACAAGATAATAATCTTTTGCTGCTTAATTTGCA
CATCCCCACCCCTTGACAACTTTAAATGCTAGTTAGGCACTTAGATGGCCCTGTTCTTGTTAAACTGC
TCTTAGCTAAGATGCAATTCTCAGTGCTTTCAAGTGGATTCTGTTGAAGAAAATCTTTGTAATAGC
CTTTTGATGCTGCTGTGTACAGTCTTCATTATGCATTGGGCAGTATTTCTGGCTAGAGTTTAAAGG
AACAGAAAGAAAACAGCTTATTTTCCTTCTTACGGACTCATCTTTAGCGTTTATTTCAACCTTTTGCT
AATTCTCTGAGAAATCTGCAGCACTCAGCCATACACCAACAGTGTTGGAAAGTTAACACCCCTGGTTAGG
GCAGAAATGTTAAAGACCATCTTGGCAGAGTTCCAGCCACGCTCTTTATTCTGTTCTCAAATAAAGCAGT
GTCAGTAGTTTTCTCTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_017880.3</u>
Locus ID:	54980
MW:	18.7