

Product datasheet for **SC206613**

DBNDD2 (NM_001048221) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: DBNDD2 (NM_001048221) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: DBNDD2
Synonyms: C20orf35; CK1BP; HSMNP1
ACCN: NM_001048221
Insert Size: 497 bp
Insert Sequence: >SC206613 3'UTR clone of NM_001048221
The sequence shown below is from the reference sequence of NM_001048221. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGAGGGGCAGAGCCTGGAGCCTGCAGCTAGCAGTGGGCCCTGCCTACAGACTGACCACGCTGGCTATT
CTCCACATGAGACCACAGGCCAGCCAGAGCCTGTCGGGAGAAGACCAGACTCTTTACTTGCAGTAGGC
ACCAGAGGTGGGAAGGATGGTGGGATTGTGTACCTTTCTAAGAATTAACCCTCTCCTGCTTTACTGCTA
ATTTTTCTGCTGCAACCCTCCACCAGTTTTGGCTTACTCCTGAGATATGATTTGCAAATGAGGAG
AGAGAAGATGAGGTTGGACAAGATGCCACTGCTTTTCTTAGCACTCTTCCTCCCTAAACCATCCCGT
AGTCTTCTAATACAGTCTCTCAGACAAGTGTCTCTAGATGGATGTGAACTCCTTAACATCAAGTAAG
GTGGAAGTCAAGCCATGCTGCCTCCTTACATCCTTTTTGGAACAGAGCACGGTATAAATAATAACTAA
TAATAATATGCCAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001048221.3](#)



[View online »](#)

Summary: May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation (in vitro).
[UniProtKB/Swiss-Prot Function]

Locus ID: 55861

MW: 17.8