

Product datasheet for **SC206616**

RLF (NM_012421) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: RLF (NM_012421) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: RLF
Synonyms: ZN-15L; ZNF292L
ACCN: NM_012421
Insert Size: 504 bp
Insert Sequence: >SC206616 3'UTR clone of NM_012421
The sequence shown below is from the reference sequence of NM_012421. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACAGATGAGCTTTGTGTAGGAAGTTCATAAGTAGCAATTTGTTTTAGTAACAGACTGGCTCCAACACT
GCAACATGGGGACATTTGCCAACTCGAACAAAGGCTGAGAAGCAGCCACACCGTTGTTTAGGGTAGAAT
AGGCTGTGTATTTACATGAATGTATAATCTATGTCAGCAGTATTGGCTGAGTCCATTAGCTCTCCAG
TTGGTTTAAATGATTGGGTTTATTTTTGTTGTTTATTAAAAAATGGAAGTACACTTGTGTTGG
TGCTAATTAATACATCAAAATATACTGGGGCTTCTTTTTCAAATTAAGTGTGCATGATTGTATATGGA
ACAAATACTAAGGTCCCAGGGTGGGAGGGCTAGGGAAAGGGATATGGAGTTCTTACTTGACTTGAATGT
GCACCTGAGGGTGCTTTGTGTAATATATTGTACACTACAGCATCTTATTTTTTTGAGTTGAGTTTCAA
TAAATTACAATTTTTACCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_012421.4](#)



[View online »](#)

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 6018

MW: 18.9