

Product datasheet for **SC206621**

SIRPG (NM_001039508) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: SIRPG (NM_001039508) Human 3' UTR Clone
Symbol: SIRPG
Synonyms: bA77C3.1; CD172g; SIRP-B2; SIRPB2; SIRPgamma
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_001039508
Insert Size: 517 bp
Insert Sequence: >SC206621 3'UTR clone of NM_001039508
 The sequence shown below is from the reference sequence of NM_001039508. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCTACGTCCCCTGGAAGCAGAAGACCTTGACTCTCCTTCTTCTCCCTGCCACGTGGGACCCCTCATC
TCTGCTGCCTCCTTCTTCTGAGAGGCTCAGCTTGAGAGAATGAGCCAGTGAGAAGCTTCTCTAGAC
TTGGCTCAAACATCTCCCCTCCAAGACATCTGCCTGCCACAGGCTCCTGTTGCTCCTTCACACAGA
CCTGGATGCCCCAGAGCAAGGTCTTCATTCATGGTCCTGAGCAGGTGCCATGGGATTGGGCTCTGGGCA
CTGACTTAACGGCACCTCCCTAGAAGGGCGAGAAACATGCCAAATCTAAACACACCAGGACTCCCATCCA
TCGCCTTGAGACTGACCGTAAACACAGACGCTCTCCAGTTCTCAAGAGTTATCCTGCCTTCCAGATT
CCTGCCTATCCCAACTCCCAGCCTTGTGAGGTTCTCTATTGCCTCTTGAATACAAATGCACTCCCAA
AGTGGTTTTAAGAAAATAAAAAGATTATCCTTCC
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001039508.2
Summary:	The protein encoded by this gene is a member of the signal-regulatory protein (SIRP) family, and also belongs to the immunoglobulin superfamily. SIRP family members are receptor-type transmembrane glycoproteins known to be involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]
Locus ID:	55423
MW:	19.1