

## Product datasheet for **SC206639**

### N acetyl transferase 5 (NAA20) (NM\_181528) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** N acetyl transferase 5 (NAA20) (NM\_181528) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** NAA20  
**Synonyms:** dj1002M8.1; NAT3; NAT3P; NAT5; NAT5P  
**ACCN:** NM\_181528  
**Insert Size:** 543 bp  
**Insert Sequence:** >SC206639 3'UTR clone of NM\_181528  
The sequence shown below is from the reference sequence of NM\_181528. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TATGAGGAAAGCACTTTCCAGGGATACTGA GAAGAAATCCATCATACCATTACCTCATCCTGTGAGGCC
TGAAGACATTGAATAACCCTGGGCAGTGGTCTTAGGCAGATACTCTAGATGCTTTATGGACAATATTA
TTTTATTGGATGATTCTGGAGCTCTATTAGGAGAAAAGTAATCATTTTTAGGTCTTAAAGACTTCAAGA
AAATACAGGTTATCAATTTATTTAAATCTCATTGTTTCCAGTTAGCAATATCATACCTATTAAGCTG
TTCATTGTAACAAAATCAATCAAAAAGGCAGCTAGGTCAGAAGGAAACATACCACTCTCATGGTTCAT
AGTATTCACTGTATGTATGCTAGGAAAAGACTTGCTCCAGTCTCCTCCTCAGTTCTGTGCCTGAGAAC
CACTGCTGCATATATTTGTTTTAAATTTGTATTGAACTGTTAATTGAAGCTTTAAAAGCATATATGA
AATGTATAAATCTAAGATGTATAATACATTATTGACTCTATGAAAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_181528.3](#)



[View online >](#)

**Summary:** NAT5 is a component of N-acetyltransferase complex B (NatB). Human NatB performs cotranslational N(alpha)-terminal acetylation of methionine residues when they are followed by asparagine (Starheim et al., 2008 [PubMed 18570629]).[supplied by OMIM, Apr 2009]

**Locus ID:** 51126

**MW:** 20.9