

Product datasheet for SC206669

C6orf115 (ABRACL) (NM_021243) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C6orf115 (ABRACL) (NM_021243) Human 3' UTR Clone
Symbol:	C6orf115
Synonyms:	C6orf115; Costars; HSPC280; PRO2013
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_021243
Insert Size:	429 bp
Insert Sequence:	>SC206669 3'UTR clone of NM_021243 The sequence shown below is from the reference sequence of NM_021243. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GATGTTGACATTATTAAGTGAAGATAATGTGGTTTACATATCTTTATGTAAGTACTGCAATTTTTGTTT
CTGGTAAACTGGAATATAAAGTGAAGAAACAACATTTGAACATACTTAATGTATTTTTATAGAATTT
GTAACCGAAAGGAGATTCATGTTTTAGAAGTCTGTCCTTTTTTATATCTTGAAAGAAAATCTATGTATG
ATGCTATAAAATAAATCCTATTATTTTTCTCAGGAATCTGGTTAGGAATTGCAGGCAATGAGATTTTTT
GCGGGGCAGGGATGGGAATGTTTGTTCATAAATAATTAGACATTTCTATAGATATTTGACATTCTGCG
AAAGCAACAAGCAAAGTGAAGACCAACTCCTATGAGAAATATTATGATGTTTATGTAATAAGACATGT
AACTGTCTTAAATTT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_021243.3](#)

Locus ID: 58527

MW: 17.1