

ZCCHC9 (NM_032280) Human 3' UTR Clone

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_032280.3</u>
Summary:	May down-regulate transcription mediated by NF-kappa-B and the serum response element. [UniProtKB/Swiss-Prot Function]
Locus ID:	84240
MW:	25.8