

## Product datasheet for **SC206706**

### CCDC83 (NM\_173556) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC83 (NM_173556) Human 3' UTR Clone
Symbol:	CCDC83
Synonyms:	CT148; HSD9; KP-CoT-23
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_173556
Insert Size:	533 bp
Insert Sequence:	>SC206706 3'UTR clone of NM_173556 The sequence shown below is from the reference sequence of NM_173556. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACATATAAGATGATGAAGTCTTTTCTCTAAGACGGAAAGCTGCAAAGGAAACACAACCTTTTCCTTATAA
ATGTTCTTTGGAACTGAAGTATATCCGTTGCCATTTTACTTACACTTTGGCTCATTTTTAAACCAGC
TGTTATTTCTAAAGTCAATTTACATTTAAAAATCAAAGGTATTCAGCTATTCATTTACTTGCATGGTA
TGAGTGACAAAACGGAAGCACGCTTTGTATTTCTACACTGAAGTATTCAGAAGCATGACAGTGGGTTTC
AAGGTAGTCTCTGAGTTCTTTTACACACAAAAAATCACTGATTAATCTGTGATTCCAGTATGAAA
TAGTTCATTAGAAATGTTTCTAAGAAAACTTAGAAGTTTGCATAGCATTGTCTACACATCTTTCCCT
CTGAGGATGCTCAATGTGATAGACAGCCAGTCTATAATGCAAGCCAATTCTCCGTAGTTTAAACCCTGTG
TATTAGTCTGTTCTCATGCTGCTAATAAAGACATAATTGAAACTGGGTAA
ACGCGTAAGCGGCCCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u>NM_173556.5</u>
<b>Locus ID:</b>	220047
<b>MW:</b>	20.7