

Product datasheet for **SC206765**

ZNF880 (NM_001145434) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF880
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001145434
Insert Size:	511 bp
Insert Sequence:	>SC206765 3'UTR clone of NM_001145434 The sequence shown below is from the reference sequence of NM_001145434. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCCATACTGGAGAGAAACCGTACAGATGAATGTGTGTGGTAAGATCTTTAGTAATAATTCACACCTT
GCACAGCATGAGATAATTCATTCATGAGAGAGTTCTTACAACTGAGTATGGCAAACCTTCATGCTAA
GTTCTAGCATTAAATCAACATCAGAGATTCCATACTAAAGAGAAATCATAGCAATGTATGTGAATCAGGT
CTCTTGAGGCCTGCCAATGACTAGATATCAAACATACATCTTGGATGAAACCATACAGATGGATTAT
GTATGCTGAGGCTATTATTCAAGGACCATTACTATGGAACATGATAAGATTTACACAAGCGATAATTCA
GTCTCTAGTTCTCCAATTTTATGATACTGCATGCTGCAGAAAAGTCACAGCTCCAAATACGTTGAATC
TCATGACGCGATGCATGTCGAAGGTCCAAGGCCATGTGGAATGGTCAGTTATATTCAGGCAATTAATA
AACCTAATTATATGGGGTTTTTTGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001145434.2</u>
Locus ID:	400713
MW:	19.6