

Product datasheet for **SC206903**

JUNB (NM_002229) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	JUNB (NM_002229) Human 3' UTR Clone
Symbol:	JUNB
Synonyms:	AP-1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_002229
Insert Size:	530 bp
Insert Sequence:	>SC206903 3'UTR clone of NM_002229 The sequence shown below is from the reference sequence of NM_002229. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTGTAACAATTG
GCAGAGCTCAGAATTCAGCGATCGCTTGGCGCGCC
CTGCTTGGGGTCAAGGGACACGCCTTCGACGTCCCCTGCCCTTTACGGACACCCCTCGCTTGGAC
GGCTGGGCACACGCCTCCACTGGGGTCCAGGGAGCAGGCGGTGGGCACCCACCCTGGACCTAGGGGC
GCCGCAAACCACACTGGACTCCGGCCCTCTACCCTGCGCCAGTCCTTCCACCTCGACGTTTACAAGC
CCCCCTTCCACTTTTTTTGTATGTTTTTTCTGCTGAAACAGACTCGATTATATTGAATATAAT
ATATTTGTGATTTAACAGGGAGGGGAAGAGGGGCGATCGCGCGGAGCTGGCCCCGCCCTGGTAC
TCAAGCCCGCGGGGACATTGGGAAGGGGACCCCGCCCCCTGCCCTCCCCTCTCTGCACCGTACTGTGG
AAAAGAAACACGCACTTAGTCTCTAAAGAGTTTATTTAAAGACGTGTTTGTGTTGTGTTGTTGTTT
TTTTATTGAATCTATTTAAGTAAAAAAAATGGTTCTTTATTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	AscI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_002229.3</u>
Summary:	Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'. [UniProtKB/Swiss-Prot Function]
Locus ID:	3726
MW:	19.8