

Product datasheet for SC206952

UPK3BL1 (NM_001114403) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: UPK3BL1

Synonyms: UPK3BL; UPLP

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_001114403

Insert Size: 530 bp

Insert Sequence: >SC206952 3'UTR clone of NM_001114403
The sequence shown below is from the reference sequence of NM_001114403. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCAGCACTCTGCCGAGGGGGCTTCCGAGGGGTTCCAGAGGGGCCACGTGTCCCTCCACCTCCTCC
CTGGCCCGAGGCTGCAGAGCCTGAGCTGGGACACGCCCTGAAGCTTCTGGACCCTGAGAGAGATTGTTTC
TTTTCTTTACTATTTTTTTTTTTTTTTTTTTTGTAGACGGAGTCTCACTCTGTGCCAGGCTGGCGCGATC
TTGGCTCACTGCAACCTCTGCCTCACGGGTTCAAGCGATTTTCCTGCTTCAGCCTCCTGAGTAGCTGGG
ACTACAGGCACACGCCACCATGCCAGCTCATTTTGTATTTTAGTAGAGATGGGGTTTACCATGTT
GGCCAGGATGGCCTCCATCTCTTGACCTTGATCCGCCCAGCTCGGCCTCCCAAAATGCTGGGATTAC
AGGCGTGAGCCATCAAGTCTGGCGAGAGAGATTGTTTCTAGATGAGGGTGGGGCGGGTGTCTTAGCC
CAAAGCTTGTGCCAGTCTCTATCGGAAATAAATGCCCCCAAAACCTC
ACGCGTAAGCGGGCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001114403.3</u>
Locus ID:	100134938
MW:	19