

Product datasheet for SC206999

Synaptotagmin 3 (SYT3) (NM_001160329) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Synaptotagmin 3 (SYT3) (NM_001160329) Human 3' UTR Clone
Symbol:	Synaptotagmin 3
Synonyms:	SytIII
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001160329
Insert Size:	538 bp
Insert Sequence:	<p>>SC206999 3'UTR clone of NM_001160329 The sequence shown below is from the reference sequence of NM_001160329. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGACTATCAGAGAAAGAGAACTCCGAGTGAGGGGTCTGGCCTAGGCCCGGGATCGGACCAGGCTCCCTC AGGACCCATCCTTTCTGCCCGACCGTGAATTCATCTCCTTGAAGCCATAACGTCGAGCTGCTGGT GCGGGGCAGCCCTGGCCCTAGGCTTCTAACCCCTGGAAGCGAGAGATGAGAGGAGGCCGCCCAGCTC CTTCTTTCAGGGTGGGGTCAATCAGCCTCCACTGTGTCTGTCTTTTCTCCCTGGGGCTCCCCCTCGA GGCGAGGGGCCATGCATGTCTGGGGACCCCTGCCCCCAAAACCCTCTGTCTGTCTGTCTCTTTGC TGTTTGTCCAAGACTCAGTGTCCCGACCCCTTGTCTCGCCGTGAATGTCAATGGGCCAATCCTCTCTGT CCTTTCAGACACACACACCTGTGTCCACCCCTTCTGTTCCGCCACCCCTGCGTCTGGCCGGTCCCC CCACTGTCTGTCTATCAACGCCAGAATAAACACACTCTGTGGGTCTCACTCAA ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001160329.2
Summary:	Ca(2+) sensor involved in Ca(2+)-dependent exocytosis of secretory vesicles through Ca(2+) and phospholipid binding to the C2 domain. Ca(2+) induces binding of the C2-domains to phospholipid membranes and to assembled SNARE-complexes; both actions contribute to triggering exocytosis (By similarity). Plays a role in dendrite formation by melanocytes (PubMed:23999003).[UniProtKB/Swiss-Prot Function]
Locus ID:	84258
MW:	19.4