

Product datasheet for SC207015

ASB3 (NM_145863) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ASB3 (NM_145863) Human 3' UTR Clone
Symbol:	ASB3
Synonyms:	ASB-3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145863
Insert Size:	553 bp
Insert Sequence:	<p>>SC207015 3'UTR clone of NM_145863</p> <p>The sequence shown below is from the reference sequence of NM_145863. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p>

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CCAGAACTGGCAGCTATTCAAGATGGAATATCAGTGAACTACTTAACACAGCTAATTTTTTCTCTG
AAAAATCATCGAGACAAAAGAGCCACAGAGTACAAGTTTTATGATTTTATAGTCAAAAGATGATTATT
GATTGTGAGATAGGTTAGGTTTTGGGGGCCAGTAGTTCAGTGAGAATGTTTATGTTTACAACACTAGCCT
TCCCAGTAAAAAAAAAAAAAAAAAAAAAAAAAATTGTAACATCACTTATATTACTTTATTGCAGCTTC
ATCACCAGTACATTATATGTTGTAATATTTATTTACCTGATCATTTTGATCATTTTCTGCTTTATTTTG
CTAATAAACTGTGATGTTACTTCTAGTGCTAAACATGGCATATTTCCACCTATGATTCGTGTTTACCTG
GTATTAGGAGCTCAGAAATGGAATGCATAAAGCTTCACTGGAAGTGATACAACCTGTGGTGTAGAATCTG
TTATTATTATCATTATTATTTTATTTAGACTTGACTATCTCTTATGTTTATTAAGAACATGTTTTCCTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_145863.3
Summary:	The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Jan 2011]
Locus ID:	51130
MW:	21.2