

Product datasheet for **SC207019**

ZNF443 (NM_005815) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF443

Synonyms: ZK1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_005815

Insert Size: 480 bp

Insert Sequence: >SC207019 3'UTR clone of NM_005815
The sequence shown below is from the reference sequence of NM_005815. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCCTTGCATAGACATAAAAGACTCACAGCATTCTCTCTAAATGTATGGAATGTGGGAAAGCATTAT
TAATTTATTTTCATTTCAGATACTTGTACGAAACACATTGGAGATAGACCCTGTGAATGTAAGCACTTG
GTAAAGCCTTAAGTAGTTTCAGTTTCTTTCCAGTACAGTCATCCCTTGATACCTGCTGGGTATTGGTTC
CAGCACTCCGTGAGCCATGTCCAGTCCCTTTTATAAAATGACATGTTTGTATGTAACCTACCCACATCC
TCTTGATAGTCTCAATTATGTCTAGATTACTTAAATACCTCATGCATTATAAAAGCTATGCAAAATAG
TTGTTCTATTGTATTGTTAGGGAATCATGATAAGGAAAAGAGTCTATGTATTTTCAGTATAGATGCAG
TAATTGTGAGCCTATCAACATAGTGTAGTCAGCCAGAACATTAAAGTTTCTTTGTTTCAACTCTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_005815.5
Summary:	Zinc finger proteins (ZNFs) bind DNA and, through this binding, regulate gene transcription. Most ZNFs contain conserved C2H2 motifs and are classified as Kruppel-type zinc fingers. For a general description of these proteins, see ZNF91 (MIM 603971).[supplied by OMIM, Jul 2002]
Locus ID:	10224
MW:	18.6