

Product datasheet for **SC207087**

ZCCHC18 (NM_001143978) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZCCHC18
Synonyms:	PNMA7B; SIZN2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001143978
Insert Size:	549 bp
Insert Sequence:	>SC207087 3'UTR clone of NM_001143978 The sequence shown below is from the reference sequence of NM_001143978. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GAACACAGTGATGCTTCTGAGCCACAGTAAGGATCTAGTCCAGCCCTAAATGAGTCCTTGACTGTATTC AGAGTCTGGTAATGGGAATAACAGGAGAGGGGGTGGGTTTCTAACTGCATGAATTAATCCACAAAGCA GTTTTCTTTGGGAAGGAGAAGAGTCTTGCATACCAGCACAGTGGATGGGGAATGGTGTGACCTCTGT CCCTGCTCTCTTCTGCTGGCTTAAGGGTCTATTCTCCATGTGTCTATTCTTTGATGACTTTATACT TTTTATGAAAGGAGCATCTTTTCATTAAACCTTCAAAAATAAAGGAAGATACAAAACTAAAGTACAAA TTACCTGAACTCAGCACCTTTTGAACCATTTGTCAGCCCTTTGGTGAATGTCCTTCAGATATCTCC CTGTACCTGTGTACCCAGATAGATATATGTATAGATAAGAGTGACCAAAATATAAGTGCTGTTCTATCT GTGTATTTTACCAAATAATGTATCTCGTGGACTTCTATGTCGATGAATATATATAAATATCTAA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001143978.3</u>
Locus ID:	644353
MW:	21.2