

Product datasheet for SC207112

ZNF700 (NM_144566) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: ZNF700

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_144566

Insert Size: 540 bp

Insert Sequence: >SC207112 3'UTR clone of NM_144566
The sequence shown below is from the reference sequence of NM_144566. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGTAAGGATTGTGGAAAGCATTACAGTACCTGGTTCCTTTTATGGACATGAATAGACTCACACTGGA
AGGAAGCACTATGAATGCAAGCAATGTGGCAAACTTTTACATTTTCCAGTTCTTTTCGATATCATGAA
AGGACTCACACTGGGGAGAAACCCTATCAATGTAAGCAGTGTGGGAAAGCCTTCATTCTTTTACTTCT
TTTCAATGTCATGAAAGGACTCACACGGGAGAGAAACCCTATGAGTGTATTCTAGTTCGGTTTGATATC
ATGAAAGGACTTACACTGGAGTGAAACCCTATGAATGTAAGCAATGTGGGAAAGCCTTCAGATGTGCCT
CGCACCTTCAACGGCATGGAAGGGTTACACTTGGGAGAACTCTATGAATGTAAGCAGTATGGGAAAG
CCTTCAGATCTGCAAGATTCTTTGAATACAGATAATTAATGTAAACAATTATCATAAGTATACTAACA
TGTTATTCTTTTAAATAAGAAGGTATAATAAAATATCCATTGGTTTTATGTATTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_144566.3</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	90592
MW:	20.8