

Product datasheet for **SC207118**

ABI3 (NM_016428) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ABI3 (NM_016428) Human 3' UTR Clone
Symbol:	ABI3
Synonyms:	NESH; SSH3BP3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016428
Insert Size:	540 bp
Insert Sequence:	>SC207118 3'UTR clone of NM_016428 The sequence shown below is from the reference sequence of NM_016428. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCTGGGAACTATGTGGAGCCAGCTGCTTGACAGCCAGGGCTCTCTGGGCAGCTGATGTCTGCACTGAG
TGGGTTTCATGAGCCCCAAGCCAAAACCAGCTCCAGTCACAGCTGGACTGGGTCTGCCACCTCTTGGG
CTGTGAGCTGTGTTCTGTCTTCTCCATCGGAGGGAGAAGGGTCTGGGAGAGAGAATTTATCCA
GAGGCCTGCTGCAGATGGGAAGAGCTGGAACCAAGAAGTTTGTCAACAGAGGACCCCTACTCCATGC
AGGACAGGGTCTCTGTGCAAGTCCAACTTTGAATAAAACAGATGATGTCTGTGACTGCCCCACAG
AGATAAGGGGCCAGGAGGATTGAAAGGCATCCAGTTCTAAGGCTGCTGCTAATTACAGCCCCAACCC
TCCAACCCACCAGCTGACCTAGAAGCAGCATCTTCCATTTCTCAGTACCCACAAAGTGCGAGCCACA
TTGGACCCAGACACCCCTCTGCAGCCATTGACTGCAACTTGTTCTTTTGCCATTG
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_016428.3</u>
Summary:	This gene encodes a member of an adaptor protein family. Members of this family encode proteins containing a homeobox homology domain, proline rich region and Src-homology 3 (SH3) domain, and are components of the Abi/WAVE complex which regulates actin polymerization. The encoded protein inhibits ectopic metastasis of tumor cells as well as cell migration. This may be accomplished through interaction with p21-activated kinase. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]
Locus ID:	51225
MW:	19.8