

Product datasheet for **SC207127**

LRRC61 (NM_001142928) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	LRRC61
Synonyms:	HSPC295
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001142928
Insert Size:	545 bp
Insert Sequence:	>SC207127 3'UTR clone of NM_001142928 The sequence shown below is from the reference sequence of NM_001142928. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCGGGCCCCACCTCTTCTTCGTCTTCTGACGTGGCCTATGGCCAGGACAGCCTGGCAGGTGGCCTC GCTGCCCCCAGTTCCTCTCTGCCCCCACTCGTCTTAGTTGCTTCACACTGGTCACTGGCCCTGCA CACTGGGCTATTGCTTTATCCCTATCCTGAGAGCAGCCCTCCCCACCATCCCTCCACATGCTGCAAGG ACAGACTGAAGGGCTGTGAGCAGGTGTAAGGGCTCCACATCCGTGAGCCTGTGTCCGAGCTGCTGCC ACTCTGGGCTGCTCCAGCCTGCAACTTAGTGGAAGGAATTACTTCTCCTGAGGCTACAGGCGAGAAAG GTAGGGATGGGCCAGCCTCCCGTCTCAGCTGTTGGGAGACAGTAGGCAGGCTGAGTGGCCAGAGCACT CCTGGAAGTGGGACTCCCTGCCTTGCAAATGTGCCTCTCCAGACTGCTCCTGCACTTACCCCTCCCC GCAGCACCTTCTCTGCCGTTCTTGTCACACATCTATTAAATGCTTCTGTTTTCAATTGCA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001142928.2</u>
Locus ID:	65999
MW:	18.8