

Product datasheet for **SC207205**

RNF214 (NM_001077239) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	RNF214
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001077239
Insert Size:	1333 bp
Insert Sequence:	>SC207205 3'UTR clone of NM_001077239 The sequence shown below is from the reference sequence of NM_001077239. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACTTGTCCCTTTTGTCCAACCTTAAATACGCGGACCTGACTGGGGAGGAAGAAGAAGAGAACTGATGT
GAACAGGAAGCGCGGGTTCAAGATTCTAAACTCTATATTATACAGTGACATATACTCATGCCATGT
ACATTTTATTATATAGGTAATGTGTATAGAAAGTCTGTATTCCAATGTTTCGTAATGAACTATGT
ATATTATGCAGAAACAGTCTGTTCCCTCATCTTGCAATTCCTTTGGGGGATGCAGATTGTAGGGAAG
ATGATGTTTAGTTTGGCCTTGAAATATGATATCCCTGCCAGGGCTGTTTCAAATACAATATAAAAA
CCACCTAGGAACCTGCTGTTGCTCTAAGGCCATTCTGCTTTGGTTTGGCTCAGCCTCTAGTCCATTTC
TTAAGGCTCATGTATGCAGATTAAAGCCTGGTGTCAACCACTGTCCAACAGATGCCTTGCTTACCG
AAAGCCTCCAGAAGCCTCAGTATTGTTTATGCCACTCTACTCCAATGGATAAAATGAGACTCTGATTG
AGGAAAAAAGTAACCTAGTAGTTTGACTTCTTATTGAATTTTACTGTGTTAACGCTCATTATTT
ATACAGACATTAGGTTTACAGAATATTCTGTTTACATCACCAAAATTCACAGTCCGAGAATAACAACA
TAACCAGGTCCCAATTCTCCATGTACCCACAAAGCTTCTGTCCACCTATTTCTGGACAGAAATTAG
CACAACCCACAGGTTTTCTGGGCCAAGTCTTCTTTGCTGCCACTGTCTGGCTTCTAATCAAGCTC
TGACAGGCCAACATTGTGAAGTCTCACCTTTCCATTCACTTCTGGTCTCTAGTCTAGCTAATCCC
CCTCCCCCAGAAGGTTAAGGCAGTTGGTTCAACAACCAAGTTGTTTGGTGACCTGGGTGAGCTCCAGAC
AGGCACTGGGTGAGGAGATGTCTGTGCAAAATTAATTGTAGAATGATCTTGGGATGCTGGATTTAAT
TTTTGAAGCAAGAAAGTAGGGGGGACATATTCCTGTCTGGGAAAAATGGAGACCAAGTGATATT
AAGAAAGTAGTTTGGATTTCAAATTTATAATCACTTGCTCTGTGGTAGCTATCAGGAAGTAGATCAT
TTTCCACCTCTGCTTATTGTACTTCTAATTTTAGGTTTATACCTGGGAAAGCTGTGACTCTATTAGA
GCTTTGAATCTTTCTATCTTTTATCTTTACTCTACAGGAATGGCCTCAGCTCCCTTCTCTTTTA
GAATTAACAAGTAAAACTAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```



Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001077239.2</u>
Locus ID:	257160
MW:	51