

Product datasheet for **SC207252**

PAMR1 (NM_001001991) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PAMR1 (NM_001001991) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PAMR1
Synonyms: DKFZP586H2123; FP938; RAMP
ACCN: NM_001001991
Insert Size: 559 bp
Insert Sequence: >SC207252 3'UTR clone of NM_001001991
The sequence shown below is from the reference sequence of NM_001001991. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AAAGACTGGATTGAAAGAAATATGAAATGAACCATGCTCATGCACTCCTTGAGAAGTGTCTGTATAT
CCGTCTGTACGTGTGTCATTGCGTGAAGCAGTGTGGCCTGAAGTGTGATTTGGCCTGTGAAGTGGCT
GTGCCAGGGCTTCTGACTTCAGGGACAAAACCTCAGTGAAGGGTGAGTAGACCTCCATTGCTGGTAGGCT
GATGCCGCGTCCACTACTAGGACAGCCAATTGGAAGATGCCAGGGCTTGAAGAAGTAAGTTTCTCAA
AGAAGACCATATACAAAACCTCTCCACTCCACTGACCTGGTGGTCTTCCCAACTTTCAGTTATACGAA
TGCCATCAGCTTGACCAGGGAAGATCTGGGCTTCATGAGGCCCTTTTGGAGGCTCTCAAGTTCTAGAGA
GCTGCCTGTGGGACAGCCAGGGCAGCAGAGCTGGGATGTGGTGCATGCCTTTGTGTACATGGCCACAG
TACAGTCTGGTCTTTTCTTCCCATCTCTTGTACACATTTTAATAAAATAAGGGTTGGCTTCTGAAC
TACAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 µg dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_001001991.3](#)

Summary: May play a role in regeneration of skeletal muscle.[UniProtKB/Swiss-Prot Function]

Locus ID: 25891

MW: 21.2