

Product datasheet for **SC207267**

SFMBT1 (NM_001005158) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SFMBT1 (NM_001005158) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SFMBT1
Synonyms:	RU1; SFMBT
ACCN:	NM_001005158
Insert Size:	1599 bp



[View online »](#)

Insert Sequence:	<p>>SC207267 3'UTR clone of NM_001005158 The sequence shown below is from the reference sequence of NM_001005158. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTTGCTTTTTATGAGCAGTTTGCCAACAGAAGGACAACCAAGTGAGCTGGATCTTTGAAGCACAAA TGCAGCAAATCCTTCACCTGCTTTATAAGTGGAGCTGGAATAGTCTGGGGCTCTGGGGCTGCAGGT ATCAGCTTGCTCTCTTTGCACTTTGGGGAAGGAGGACTCACAGTGAGGAAGCAAAAAGTGCACAGA AGTGGATCACCTGCTGGTGGAAATGTGGACATCTCTTGTTCAGCAGATGGCAGTTTTTAAAAATAAAG GTTGTGAGGAAAAGACTTATAAAGAAGAAAAGCATTTCAGTGGTGTGGCCTGAAAACAAAGAATAAC CTAGGCTGCTGGAAGCACCCCTTTGGTGTGTTTCATTCTGTTCCCTCCCATTGTAGATTGAACCTTGT TCTCTGCTTTCTTTTCTGGAAAGAGAGGACTTAGCTTAAAGTCAGCACTGATTTGGGACTGTTCTTA AGGCATATCAGTCTTCATTGTGCTTTTAAACTTTTTAAATTAACAGTTTCATTTGGGGA TGATTATGTGGCTCTAGGGTTTTGAATGTATAGTCACACTTTGATCATTTTAAAAATCGATTCTTAGGAG TTCTTTGTTTACTACCTCTGTTGATAATTGAGCTTACAGCCATGTATGAGGTTTTTTGTATGATGCCA TTTTTAAGCTACATAAACTAAAATTATGACAGAATTTGGAGGGGGGAAACTGTTTGCTTTCTAAA AACATTATATGGGGCACTTACACCAGTTCCTTAACTGACCTGATCTGAGAAGTCTACTGTTACTGTTCT CCATCTTGACCCAGCCGAGGCTGCCTGTGTTAAGTAACATCACCAAGACTGTTAGTATTTCTGTTAAT TGCCCTCTTCTATTCTCAGAGGGACTTCAGTGACTTTTATTTGTGTTTACTTTTCTCTAGAGT TGTTGTCTACCTCTGTTTCATTTCTAAGGCTTTACTTTGGGTGAAAATCAGTAAGTGTCCATGGAC CCACAGCGAAGCCCATTTGCTCTTTCCAGGTGTCTCCAGTGTAGAAGCACAGCACAGAAAACATGT TTTGAGGGTAACGCAGAGCATGGAGTTTCTGTTAGCAGTGTCCCTCAGTGGCAGGATAAACTAGAGAT GTGGTTAAGGATTCTGGAACCTCTTTGCTAGCCACACTTAGACTTTTTTCAGTAATCCTTTAAATT GAATATAATCTTTAGTCCGAGGATTCCTTTAAACAAGTATTCTTTAAGTGTGCTTTTAGCC CAGCTCAAGTTAACATGAACTCCAGGTGAATTTCACTTACTATTCCATTGCTAGCTATTATAGAATA ATAGCGTTTGCTGCACTATTCTAGTCATGAGGAGTGAGGAAAAACATACATAGGTTGTCTCCCTAAAT GGTTTGCTTTGCTGAATTTTTTTATCCTAATCACCTGCAAAATTATAACCAATCATTCTCTCAAG ATTTGAAAATGTTTAACTCTGCAAACTGGGTGAATTCTAAACAGCTCTACAAATATAAGGAAATATA CAGCAAAATTTTA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001005158.2
Summary:	This gene shares high similarity with the Drosophila Scm (sex comb on midleg) gene. It encodes a protein which contains four malignant brain tumor repeat (mbt) domains and may be involved in antigen recognition. [provided by RefSeq, Jun 2012]
Locus ID:	51460

MW: 61.2