

Product datasheet for **SC207270**

VSIG4 (NM_007268) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	VSIG4 (NM_007268) Human 3' UTR Clone
Symbol:	VSIG4
Synonyms:	CRlg; Z39IG
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_007268
Insert Size:	555 bp
Insert Sequence:	>SC207270 3'UTR clone of NM_007268 The sequence shown below is from the reference sequence of NM_007268. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTGGCCACTGAGGGCAAAAGTGTCTGTAAAAATGCCCCATTAGGCCAGGATCTGCTGACATAATTGCC
TAGTCAGTCCTTGCCCTTCTGCATGGCCTTCTCCCTGCTACCTCTCTTCTGGATAGCCCAAAGTGTCC
GCCTACCAACACTGGAGCCGCTGGGAGTCACTGGCTTTGCCCTGGAATTTGCCAGATGCATCTCAAGTA
AGCCAGCTGCTGGATTTGGCTCTGGGCCCTTCTAGTATCTCTGCCGGGGCTTCTGGTACTCCTCTCTA
AATACCAGAGGGAAGATGCCCATAGCACTAGGACTTGGTCATCATGCCTACAGACACTATTCAACTTTG
GCATCTTGCCACCAGAAGACCCGAGGGAGGCTCAGCTCTGCCAGCTCAGAGGACCAGCTATATCCAGGA
TCATTTCTCTTTCTCAGGGCCAGACAGCTTTAATTGAAATTGTTATTTACAGGCCAGGGTTCAGTT
CTGCTCTCCACTATAAGTCTAATGTTCTGACTCTCTCTGGTGCTCAATAAATATCTAATCATAACAG
CAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_007268.3</u>
Summary:	This gene encodes a v-set and immunoglobulin-domain containing protein that is structurally related to the B7 family of immune regulatory proteins. The encoded protein may be a negative regulator of T-cell responses. This protein is also a receptor for the complement component 3 fragments C3b and iC3b. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]
Locus ID:	11326
MW:	20.2