

Product datasheet for **SC207302**

SMIM11A (NM_058182) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SMIM11A (NM_058182) Human 3' UTR Clone
Symbol:	SMIM11A
Synonyms:	C21orf51; FAM165B; SMIM11; SMIM11B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_058182
Insert Size:	555 bp
Insert Sequence:	>SC207302 3'UTR clone of NM_058182 The sequence shown below is from the reference sequence of NM_058182. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGAAGCAATCAGAGAAAAAGATAACTGAAGATTCTGCCATGTAATGTCAGCTTGAGTGGACTCCAG
CTGAGAAGAAAGAGAAGAAAGACTTAATTATTGAATAATTTGTCAGAGGATAAACTCCCAACCTAGACC
TTTCACTTAAAAAGTGTGAATTTGTATATGTTTTTAAAAGAACCAGTACTGGCCGGGTATGCTGGCTT
TTACCTGAAATCCCAGCACTTTGGGAGGCCGAGGCGAGTGGATCGCCTGAGATCGGGAGTTTGAGACCA
GCCTGGCCAACATGGTAAATCCTGTCTCTACTAAAAATATAAAATAGCCAGGTGTAGTGGCGCGCG
TCTGTAATCCCAGCTACTCGGGAGGCTGAGGCAGGAGAATTGCTTGAATCCGGGAGTGGAGTTGCAG
TGAGCCTAGGTGCTGCCACTGCCCTCCAGCCTGGGTGACAGAGCGACTGCGTCTCCAAAAAAAAAAGG
TAAATTAATAAATAAAAAAAAAATAATAAACCAGTATTTTGTACTAAAATAAATGCCTTTGTAAAA
AAA
ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_058182.5](#)

Locus ID: 54065

MW: 21.1