

Product datasheet for **SC207334**

TMEM69 (NM_016486) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: TMEM69

Synonyms: Clorf154

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_016486

Insert Size: 561 bp

Insert Sequence: >SC207334 3'UTR clone of NM_016486
The sequence shown below is from the reference sequence of NM_016486. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAAGGACATAAGAGACCTGGTCAAGTATAAAAATATAAAAGTCTGGGAAGTGAGGAGCACCTCTGCCC
AGCTGCTGCCCGTCTGGGAAGTGAGGAGCGCCTCTGCCTGGCCGCTGACCATCTGGGAAGTGTGACA
AGCGCCTCTGCCCGGCCGCTGTGCAACCTTCCACGTGTGAAGTGACAGCCTTGTGTGATCTTTCTG
TCTTCCCCAAGTTTGCATTTTCGACATTAAGTTTACTTTTTAGTTAAAAGTTTAAAAATATATATAT
ATAAATACACTGTAGATAACATTTGTATGCCAGCTACACCTTTTCTACTTCTGTTTGGCTTTTTTCC
CCACACCAATGGTAATTTATCTTCACAGATTGTTCTTCATTCTAGAAATTGTTACTTCATGGTAATTA
CTTGAGCAAAAGCTTGAAAATCCCTGACAAGTACTTTTCATCTCATAGTATATTAGTTTCACTCAGTC
ATTTTATGAATAATATAGTTATCCACTTAAACATTTCAATATTTAACCATCTTGAAAATTAAGATTA
AAAATCCCC
ACGCGTAAGCGGCCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_016486.4</u>
Locus ID:	51249
MW:	21.5