

Product datasheet for **SC207360**

EXOSC3 (NM_001002269) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: EXOSC3 (NM_001002269) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: EXOSC3
Synonyms: bA3J10.7; CGI-102; hRrp-40; p10; PCH1B; RRP40; Rrp40p
ACCN: NM_001002269
Insert Size: 1192 bp
Insert Sequence: >SC207360 3'UTR clone of NM_001002269
 The sequence shown below is from the reference sequence of NM_001002269. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AATGTGCAGGCTATTAGCTCCAGATTGTGAATCATAACAGGAAGTGGGAAAACCTATCCACTGGAGAT
AGTATTTGGAATGAATGGAAGAATATGGGTTAAGGCAAAAACCATCCAGCAGACTTTAATTTTGGCAA
CATTTTAGAAGCTTGTGAACACATGACGTGAGTCAAAGAAAACAGATCTTCTCCAGATTGGCAGAAAG
TTGATATAGGTGGACTTTTTTACAGGTCAGTTGAGGCAAAAACCTATGGGTTTTTTCAGGTGAACCTCC
CCCATTTAAATACTCAGAAGATAAGGTGTGAATGTATGTATTATTAGAGTCCGAAAGTATTTTTATAAG
TACTGGTTTTACCCACGCTTTTGTGGGAGAGAAAATCATTGCAAAATCATTTTTTTTGTTTCGGTACA
ATAAAGTTTACTAAAAACAGTATTAGGGTTTATTTGGAGTAACTCCCTTAGTAATTGCTACAGCTTC
CATAATTAAGTGGAGTCCAGCAGTAAAATAATTTGTGATTATAACAGTAAGCCAGGTTCCACCAGAGGA
AATGTGGAATACCTAAGTATAAAGCAGCAAAGAAATGTCACCTGTAGTCCATTATTCAGCGGTGAGT
ATATAACAACGAAATACTTTTTCAATTTACCAAGTATTGTGTTTTGCATGGATTAACCTCTTTATTCTAA
CTGCAGAAGACACTTATCCTCATAATTGAACTGAGGCACAGAATTTAACTTGGCCTACATTCTCTATA
TTCATCAAATTGAAATTATAGTTATTATAGCTTATATGCCACTTTGTAATATACTTTTGTCTACTTCAA
ACAGCATGATACATTTCAACTTTAGTCACATATTTTCATGTTACAGGATATTTTGGAGTTTATCATT
ATATAAATAACACTATGAGGATTTCTGAATCCAAATCTGGCCAAGTCTTATTTTGGCTAAAATCATTTC
TAGGAGCTAAGCATAGAAAACCTTCTATGCTATTTAGTATATTAACCAAAACCTTTAGAGGTAAAGAG
ACTGCCAGGTGATCCACATCAGTATTCTTGGTTTCTAGACCTCTGAGAATATGTGAGATACAGATCCC
CTGAAGCCTATCCATAGAAAGTCCTTTTTGGATCAAAGAGATTATAGAATCTATTTGTTGATAAAGAC
CAAAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001002269.2</u>
Summary:	This gene encodes a non-catalytic component of the human exosome, a complex with 3'-5' exoribonuclease activity that plays a role in numerous RNA processing and degradation activities. Related pseudogenes of this gene are found on chromosome 19 and 21. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2012]
Locus ID:	51010
MW:	46