

Product datasheet for **SC207382**

Angiopoietin like 4 (ANGPTL4) (NM_001039667) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Angiopoietin like 4 (ANGPTL4) (NM_001039667) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ANGPTL4
Synonyms: ARP4; FIAF; HARP; HFARP; NL2; PGAR; pp1158; TGQTL; UNQ171
ACCN: NM_001039667
Insert Size: 514 bp
Insert Sequence: >SC207382 3'UTR clone of NM_001039667
The sequence shown below is from the reference sequence of NM_001039667. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGCCCATGGCAGCAGAGGCAGCCTCCTAGCGTCCTGGCTGGCCCTGGTCCCAGGCCACGAAAGACGG
TGACTCTTGGCTCTGCCGAGGATGTGGCCGTTCCCTGCCTGGGCAGGGGCTCCAAGGAGGGCCATCT
GGAAACTTGTGGACAGAGAAGAAGACCAGACTGGAGAAGCCCCCTTCTGAGTGCAGGGGGGCTGCAT
GCGTTGCCTCCTGAGATCGAGGCTGCAGGATATGCTCAGACTCTAGAGGCGTGGACCAAGGGGCATGGA
GCTTCACTCCTTGTGGCCAGGGAGTTGGGGACTCAGAGGGACCACTTGGGGCCAGCCAGACTGGCCTC
AATGGCGGACTCAGTCACATTGACTGACGGGGACCAGGGCTTGTGTGGTTCGAGAGCGCCCTCATGGTG
CTGGTGTGTGTGTAGTCCCCTGGGGACACAAGCAGGCCAATGGTATCTGGGCGGAGCTCACA
GAGTTCTGGAATAAAAGCAACCTCAGAACA
ACGCGTAAGCGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001039667.3](#)



[View online >](#)

Summary:

This gene encodes a glycosylated, secreted protein containing a C-terminal fibrinogen domain. The encoded protein is induced by peroxisome proliferation activators and functions as a serum hormone that regulates glucose homeostasis, lipid metabolism, and insulin sensitivity. This protein can also act as an apoptosis survival factor for vascular endothelial cells and can prevent metastasis by inhibiting vascular growth and tumor cell invasion. The C-terminal domain may be proteolytically-cleaved from the full-length secreted protein. Decreased expression of this gene has been associated with type 2 diabetes. Alternative splicing results in multiple transcript variants. This gene was previously referred to as ANGPTL2 but has been renamed ANGPTL4. [provided by RefSeq, Sep 2013]

Locus ID:

51129

MW:

17.7