

Product datasheet for **SC207416**

SH3GL3 (NM_003027) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: SH3GL3 (NM_003027) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: SH3GL3
Synonyms: CNSA3; EEN-B2; HsT19371; SH3D2C; SH3P13
ACCN: NM_003027
Insert Size: 486 bp
Insert Sequence: >SC207416 3'UTR clone of NM_003027
The sequence shown below is from the reference sequence of NM_003027. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTGGAAGTGATCGTGCCTTTACCTCAGTAAATGTGTAACACAAACTCTGGACATACTTTTCGTAAGTAA
ATGAATTCACACCAGTGTGCTCTCAGTGGGTTCTGTGACATCCTTTGCTCTGACCAACTTAATG
ACTTTTGTATGTGTGCTCTCTTTATAATGTATTTTATATCACTTTAATTTGTATAAATGATTTTCTGT
CCTTGCTACATGAAAATATTTCTTTTCTTCTCCTGTCCATAAAAGTCATTGGTTAAATGATTTGCTT
CCTGTGGCTAAAATAAGTCTCACCCATTGCAGTTATGTCAACGAATGGCCTATATTCTCAGCTGCAA
TGAAATGGTAACATTTGAACTAAGAAATGCTAAATATTTGTTTCTCGACATTCCTGATGACGTCTGG
TCTTTTCTTTTATTGATTTTAAGCTTACCTGTGAATAGCCCAATAACATGACACACTGTGTTGGCA
AAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_003027.5](#)



[View online »](#)

Summary: Implicated in endocytosis. May recruit other proteins to membranes with high curvature (By similarity).[UniProtKB/Swiss-Prot Function]

Locus ID: 6457

MW: 18.5