

Product datasheet for **SC207456**

iASPP (PPP1R13L) (NM_006663) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	iASPP (PPP1R13L) (NM_006663) Human 3' UTR Clone
Symbol:	iASPP
Synonyms:	IASPP; NKIP1; RAI; RAI4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006663
Insert Size:	582 bp
Insert Sequence:	>SC207456 3'UTR clone of NM_006663 The sequence shown below is from the reference sequence of NM_006663. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGGGTGAAGCCTCAAAGGAGTAAAGTCTAGCAGGATAGAAGGAGTTTCTGAGGCTGACAGAAACAAGC
ATTCTGCCTTCCCTCCAGACCTCTCCCTCTGTTTTTGTGCTTTATCTGCACCCCTCACCTGCTG
GTGGTGGTCTTGCCACCGTTCTCTGTTCTCCTGGAAGTCCAGGGAAGAAGGAGGGCCCCAGCCTTAA
ATTTAGTAATCTGCCTTAGCCTTGGGAGGTCTGGGAAGGGCTGGAAATCACTGGGGACAGGAAACCACT
TCCTTTTGCCAAATCAGATCCCGTCCAAAGTGCCTCCCATGCCTACCACCATCATCACATCCCCAGCA
AGCCAGCCACCTGCCAGCCGGGCTGGGATGGGCCACCACCACTGGATATTCTGGGAGTCACTGC
TGACACCATCTCTCCAGCAGTCTGGGGTCTGGGTGGGAAACATTGGTCTCTACCAGGATCCCTGCCC
CACCTCTCCCAATTAAGTGCCTTCACACAGCACTGGTTTAATGTTTATAAACAAAATAGAGAACTTT
CCTTATAAATAAAAGTAGTTTGCACAGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006663.4</u>
Summary:	IASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53; MIM 191170), whereas ASPP1 (MIM 606455) and ASPP2 (MIM 602143) are activators of p53.[supplied by OMIM, Mar 2008]
Locus ID:	10848
MW:	21.3