

Product datasheet for **SC207469**

SEZ6L2 (NM_001114099) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SEZ6L2 (NM_001114099) Human 3' UTR Clone
Symbol:	SEZ6L2
Synonyms:	BSRPA; PSK-1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001114099
Insert Size:	570 bp
Insert Sequence:	>SC207469 3'UTR clone of NM_001114099 The sequence shown below is from the reference sequence of NM_001114099. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GATACGCGGGAGTATGAAGTTTCCATCTGAACCCCAAGACTACAGCTGCAGGACCCAGGACGCCCTCC
CCTCCTCATTGGGCAGAGGGAAATACGGGACCCGGTCTCTGCCTCCTGGCTGCCCTCCTCCCTGGCTG
TGTAATAGTCTCCCTATCCCACGAGGGGGCTTTGATGGCCCTGGAGATCCTACAGTAAATAAACAGC
ATCCTGCCGCCAAAGCCGCCTTCTCTAGTTGCCAAACGAGGGGCCTGCCCCCGCCCTACCGGCTTT
TGGATTCTGGGAGGGAACTCTGCCTCCCTGCAAATCTTGAGCCCTCCTGCCAGGGCACCCCTCAA
GGACTGCCCCGATAGCTCTACTGTTCCCTTGGCCACGAAGGTGCCCCCTCCCAGATGCCTGGCCCT
AGGCCTGACTCCGGCCAGGAGGGTCAAGAAGGACAAAGGGGAGAGCTGGGACAAGGCCTTGGCCCT
TCTGCCATCTCCCAACCCACAGTCTCTCCACCTTTGCTTCTGAATTCTTGTTTTTGAGCAATAAACA
GAAAATCGCCACTTGTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001114099.3
Summary:	This gene encodes a seizure-related protein that is localized on the cell surface. The gene is located in a region of chromosome 16p11.2 that is thought to contain candidate genes for autism spectrum disorders (ASD), though there is no evidence directly implicating this gene in ASD. Increased expression of this gene has been found in lung cancers, and the protein is therefore considered to be a novel prognostic marker for lung cancer. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Aug 2011]
Locus ID:	26470
MW:	20.8