

Product datasheet for **SC207540**

C7orf26 (NM_024067) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: C7orf26 (NM_024067) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: C7orf26
ACCN: NM_024067
Insert Size: 595 bp
Insert Sequence: >SC207540 3'UTR clone of NM_024067

The sequence shown below is from the reference sequence of NM_024067. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTGACCTTCCCTTCAGGCCCATCCGCTAGGCTGGCCCGTGTGTGCCTTCTGCGCTCTCGTGGACGAA
GCCTTTCGAGATGGAAGGGGTGGCCGACTCCCAGAAGAGAACCTCGGGGAAGGGGTCGGGCAGCCCT
CCCGCCGGCAGAACCGTCTTGGTGTACGGAGTCCAGGTGCTTCCACCCGGTTCGATTCTTTGACAT
GCAGATTGGATGGTGGAGGGAAGAGTCCAGCCTCTGCCGAGAGCCTGCTGCGTGCATTTTTAAAAGATG
CCGATCCTGGGAGCCTCTGTTCTCTGCGCATTTTACACACAGCCTGTGTGGCAGGAGTGTGACGGCAG
GAGCCACGGGTGCAAGCCCGTGTGTCTGGCCTCTTTCCTCGTGAAGACGATGTGTCCCGCCAGAAAAA
GTGGGCTCCTTCTGCAGCCCCGTGAGCTGAGCCCAGGCTGCGTAGTGACCACAAGCTTATGTGCAGCAC
TGCTCAGGGAGGCTGTGAGGAATCCCCCACCTCGGAAAGGAACTTCTCAGTTTTATTGGGGGTGTCT
AAATTTCTTTTCATATGTTCAAATAAATTTTTCTAAACAGTCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 µg dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_024067.4](#)



[View online »](#)

Locus ID: 79034

MW: 21.9