

Product datasheet for **SC207580**

FAM174A (NM_198507) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	FAM174A
Synonyms:	HGS_RE408; NS5ATP6; TMEM157; UNQ1912
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_198507
Insert Size:	587 bp

Insert Sequence: >SC207580 3'UTR clone of NM_198507
The sequence shown below is from the reference sequence of NM_198507. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TTGTTTGTATGCCAATCATCCTCGAAGATAGAAATGTGCCTTTTGTGAAAGAACTTTATCTTTCTACAA
TGAAGAGTGGAAATTTCTATGTTTAAGGAATAAGAAGCCACTATATCAATGTTGGGGGGTATTTAAGTT
ACATATATTTTAAACACCTTTAATTTGCTGTTGCAATAAATACCGTATCCTTTTATTATATCTTTATAT
GTATAGAAGTACTCTGTTAATGGGCTCAGAGATGTTGGGGATAAAGTATACTGTAATAATTTATCTGTT
TGAAAATTACTATAAAACGGTGTCTGATCGGTTTTTGTCTCTGCTTACCATATGATTGTAATTTG
TTTTATGTATTAATCAGTTAATGCTAATTATTTTGTGATGTCATATGTTAAAGAGCTATAAATCCA
ACAACCACTGGTGTGTAATAAATTTAAATTTCTTTACTGAAAGGTATTTCCCATTTTGTGGGG
AAAAGAAGCCAAATTTATTACTTTGTGTTGGGTTTTTAAATATTAAGAAATGTCTAAGTTATTGTTT
GCAAAACAATAAATATGATTTTAAATTTCTCTTAAA
ACGCGTAAGCGGCCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_198507.3
Locus ID:	345757
MW:	23.1