

Product datasheet for **SC207606**

UBE2CBP (UBE3D) (NM_198920) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	UBE2CBP
Synonyms:	C6orf157; H10BH; UBE2CBP; YJR141W
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_198920
Insert Size:	662 bp
Insert Sequence:	>SC207606 3'UTR clone of NM_198920 The sequence shown below is from the reference sequence of NM_198920. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC TCCTTTTCAGGTGGCCTTTTGAAGATGTAACTATTAGAGCTGGGCAATTCTCCTCCCCGGCAGACAGCT CTCAAGCAGGAGCACACCGAGGCCAGCAGTCAGCTCAGGAGGCATCAATTACAAAGACAAGAACCAGAG CTATGGGATAAGAGCATTGCATTGAAGTTTAAAGTTTACTTGTTCATGTATTCTCTAAGAGAAGAGT TAACCTTTAGATCTTTGAAAACCCAGGACTGTGAAGAAAGTCATAGTAGCAGTTGAGGCAAGTGTAGAA GTTACATCTCTATTTCCACAGAGAGGAATATTCACATATGTGCAGTGTGAGCACCACCTGGAAGGTTA TTACCTAATGGAGATCCATCTGCTGTGTGGAATCCAGTGAGGAATTTTATTGGTGGAACTGAATGTT AACTGCTGTAATCCTAGGAAATAAACTTTTAAAGAGTTTCTCAGGATGATTTAAACAAATTTTAAAA AAGAAAAGTAGAAAACATTTTATCTTTATTCAGAGGAGATGCTCAATAGATGTTTGTGAATGAATGA TTATCAAATTAAGGATTTCTATTTATGCAAATCTGCATGTTTCTTTTGTATTTCCTTTACCATACA TTATTGCAATTAATAATATACCTCACTAGCATCCCGTTA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_198920.3</u>
Summary:	E3 ubiquitin-protein ligase which accepts ubiquitin from specific E2 ubiquitin-conjugating enzymes, and transfers it to substrates, generally promoting their degradation by the proteasome.[UniProtKB/Swiss-Prot Function]
Locus ID:	90025
MW:	25.4