

Product datasheet for **SC207613**

LRRN3 (NM_018334) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	LRRN3
Synonyms:	FIGLER5; NLRR-3; NLRR3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_018334
Insert Size:	585 bp

Insert Sequence: >SC207613 3'UTR clone of NM_018334
The sequence shown below is from the reference sequence of NM_018334. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GTTATAGGTTTACCAACAAATATGTCCAAAAACCACCAAGGAAACCTACTCCAAAAATGAACAAAAA
AAAAAAGCGAAAGACTGCAGTTGTGCTAAAAACAAACAAACAAACAAACAAACAAAGTAAAA
AAGATTACTTTGAGAGAGAAGTTTAAGCTTACCAATGCTGCTCCTGACCAATGGAAATATGTACAAC
TTCAGCATTTTAAGTAAGTGGCTTCAAGGGTACTGTGGCAACCAATAAAATACTCCATTTTCTAAA
ACTTTTCATGTAACCTTTATGTCTGGACTACAGTTCAAGTGGACAAAAACATTTCTGTATTTTTTTAAG
TAAATAAGAGTAGTTGAAGTGAAGCAATACCTCCTCTGTGTGATTACACATATTAGCCACGAGTTTT
TGCAGTGACCAGATAAACTTGAATTGACACGTGGTGAATAAAATGGACAAATTCTGTAGAGTAGACAC
AGTGAGTATGTGGACCTCTTTTATAAGGAAAAATACATTTTGGATTAAATCAATTGCTTCTGTCTTGT
TTTGTTTCTAAATAAGAATAATTTCTGGGAAA
ACGCGTAAGCGGCCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_018334.5</u>
Locus ID:	54674
MW:	22.5