

Product datasheet for **SC207809**

ZNF410 (NM_021188) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF410 (NM_021188) Human 3' UTR Clone
Symbol:	ZNF410
Synonyms:	APA-1; APA1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_021188
Insert Size:	1010 bp
Insert Sequence:	>SC207809 3'UTR clone of NM_021188 The sequence shown below is from the reference sequence of NM_021188. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAAGGAGATTTAACTGAAAGACGGACATGAGCGTGGGTGCTGACTCCTGGAAGAGCAACTCTATCTGAT
CTCAAAATGCGTATACTGGGAACAGGATGCCTTAGCCACAACAGAACAGAATGAATCTTTGAAGGCA
CAAGACTCTGCTTTTGCCACTCTTCCTCTTTCCTGGTATAGAAGATGGATGTAGGAGAGCTTCTTTTCT
AACTACCATCTGATCAGACAAGGAATGAAGCAATGACTGTGGGCTGGGAACTGTACCTACCTCTCTTC
CCACTGCAAAATTTCTGGGATAGACCAAAAGTGAAATTTGATTATGTGTTGGCTGAAGTTCTTCATTCTGA
CTGTTGAGGGGAGGTTTTCTTTGAAGAGTTTTTCATCCCAGACTCAGCTGTCTTTTCACATGGATGAAA
TAATTCCTGCTACCAACAACAGAGCTTCACCAGGAAGTTGAGTTTTCAAGATGCCTTGTTGCTTTGAAG
AAGGGAGTGATGTCAATTCTCTTGTTACATTCTCCCTTTAGCAACCTGAGTAAGAGACTCTCTGCCACT
GGGCTGCAAAAAAATAAATTACTTGAATCTCCCTTGCCCGAGGCTGAGGTACTATCTTGTCCTATACA
CTTCTACTTGGTCACTTTGCTTTCTTCGTTAATGGAACATACAGTAGCATGTTAAGAGGGATTTTCATGT
TTTTGTTTTTTTAAAGTTAAAAATTACATGATGCAGAGATTACAGTTTTCTTTAAATAAACAACAG
CCAGTCAGTTCTTTGGCTCTTGTTTTATACAAAATTTGTTTTCTTAGAGATTAAGAATTTAATCTTTTC
CCTTTAAATAGTGATTGATTTACCCTTAGGATTCCATACCAAGTAAACTGAACCGAGGAGTCTTAGG
AAACAACAAGAACATCTTCATTTCTTAAGCCCAGGTGATAGTTACTCTGTCAACCAAAAAAAGACGCA
TCTGAGAAAAATGGCAATAAAAAACAGATACTTCTGAATTTTTCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_021188.3</u>
Summary:	Transcription factor that activates transcription of matrix-remodeling genes such as MMP1 during fibroblast senescence.[UniProtKB/Swiss-Prot Function]
Locus ID:	57862
MW:	38.5