

Product datasheet for **SC207823**

DNAJC9 (NM_015190) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DNAJC9 (NM_015190) Human 3' UTR Clone
Symbol:	DNAJC9
Synonyms:	HDJC9; JDD1; SB73
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015190
Insert Size:	1476 bp



[View online »](#)

Insert Sequence:	>SC207823 3'UTR clone of NM_015190 The sequence shown below is from the reference sequence of NM_015190. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCCGAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AAATCTGCTCTCAAGAAAGAAAAGAAATAATTGGAATTTTTCTCTTCAAAGGTCCTTAGGTGTAAATTGA TGCCATCGTAGGCAAGGTGCAGGCAGGATTTGAAGGCAAAAGTCAATTCAGCTCTTGAGAAAAGGTGTC TTTCCAGCCTGAATTTTTCAGATTGACTAGACCAAGCAGAATCTCTCAACCTGATCTTAGTATTTCCTA GAAAGCACTTGACATTGTGTGAGGTCTCACCTGAAGGAACTTGGTGGTGACATTTGGGAGGGTGGAGGG AGGCAGTGTCTCTTCTGACAGCACTTGCCTCCATGGATCTTCTGTACACAGAACTCTTATCTAGGATGT GGTTCTGTTTCATGCTGCTTTCTGCGATGTGCGTGTCTGTTAGAATAGGCTCTCTACCCAGCTAGAACAC CTTCCAGACACTTGCTGGACAGCTATCTTCCACATACTTCCCAGTTTACATTTGGTCTTAATGATCTTG AATAGATCCTCTCTTCATTTTACTCAGCCAGGTTTTGTACTGATGTACAGGTGTTAAATTACTTCAAGC ATTTTTGTAAGAGGTGTATATAATTCAATAAAAAAGGTAAAACATGATGATTAAGTTCTGGGGGCTTTG TAAATGATCCCACTAAAATGTGACCTAGGAAAAATATGAATGGTGTTAGGATGAGAGAAAAGGGGAAA ACAATAGCCCTGGTCAGCTTTATAATAGAGAGCCTGGTTTCCCTAGCATGAAGAGATGTATGTTGTAGT CCTGCCACTGATTACTGAAGTCCTGCCATTAATTGCTGAATGTCTTCAGTTGGGCCACTGAGCTTGTCT GAATCTGTTCCCTTTTATAAAAGAGTTACTACATCAAAACAAAGAGTGAAATCCAATTTGTCAAACCTG TATGTATTAACATGTCCAGTTTTTGGTTAAAAAAAATTGCACTGGCTTTGAGGGAAAACACACAGG GTGAGGGAATTGGGCTAAATGACTTCTTACAGGCCCTTTCTGATTCTTTAACTTTGAAAGGCAAGCCA TATTGATCCAGTTGTTATAGTGAACATCATGGTAATGGTTTGTGAGAACAATAGAGATTTTCATTTCTAT GTAGATGAGTTGGTATGAGAATATATGGAATTTTAAAGGACTGTTTAAATCTTTGATTGTAGACTAT TAAATATACCGTATGCATAAAGTAAGCCTTTAGCTCTAAGGTAAAGACGACACGTTTTTCGGTTTGTGAC TACAAATAGGTTAAAAATAGATTTTAAATTTTAAAAATATAATTTAATGCAGGTTGTTTGAAGCATC TGTCTTCATATGATGGCATTAGAACACCTTGGTATAATAAAAAGTTACCGTAATTTATGATTATTTGAA TTTATCCATTCTGAAAATTAATAAGATCTAAAACCTGGCATGACAATCAAGATTTGTATTTAGTGAAATT TAAATAAATGTAAGCCATAGTTAAAA ACCGGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_015190.5
Summary:	May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8. [UniProtKB/Swiss-Prot Function]
Locus ID:	23234
MW:	56.9