

Product datasheet for SC207826

Intersectin 2 (ITSN2) (NM_147152) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Intersectin 2 (ITSN2) (NM_147152) Human 3' UTR Clone
Symbol:	Intersectin 2
Synonyms:	PRO2015; SH3D1B; SH3P18; SWA; SWAP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_147152
Insert Size:	605 bp
Insert Sequence:	>SC207826 3'UTR clone of NM_147152 The sequence shown below is from the reference sequence of NM_147152. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGCTCTCGGGTATCTGCTGTCTCTCATAGAGATGGTGGGCATCAGACTCAGGGCTGCCTCCACGCAG
AGTCAAAGCAAGGCATCACTTTTGATGTGTGAATTCACAAATAGTGACGGAAGCTCTCACATCCTCCAA
ATGCTGTCTCTGCCTGCCGATAATGCTTGAGATTGAAAGTCTCTAATGAGCTCTTCCCCCAGATGAG
GTCACCTCAGAGTGAAGCGGGGAGAACAGAGGGCATTTCGCATGGCTTCGTTGGATGTGGCAGGAGCCCA
TCAAGGAAGGACGGGATAGAGTGGATGGGAAGGGCCTATGGCAGACCGTAGCTTCCTTGGATATTGTC
CTAATATCTGTTTAAACATCTGACATTTTCATAAACTGGTATCTCTGGAGAACTGTGAAACAGTGAAA
GTGTTACCTCATGGTTGTATCAGTTTGGAAAACCAATGGGAGCATATTGTAAATAGTTCCCAAATA
TCATATAGCTACTGTTTGTATACCAAGTGACCTCTACGCTGATGACAGCTATCCTTATTCGAGTAGCA
CTTAAATGATTGTGCTTGAGTGAACAAAAGAAGACTTTCCATTTCTACTA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_147152.3</u>
Summary:	This gene encodes a cytoplasmic protein which contains SH3 domains. This protein is a member of a family of proteins involved in clathrin-mediated endocytosis. Intersectin 2 is thought to regulate the formation of clathrin-coated vesicles and also may function in the induction of T cell antigen receptor (TCR) endocytosis. [provided by RefSeq, Jan 2017]
Locus ID:	50618
MW:	22.3