

Product datasheet for SC207842

C1orf84 (NM_001012961) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C1orf84 (NM_001012961) Human 3' UTR Clone
Symbol:	C1orf84
Synonyms:	FLJ34502; RP11-506B15.1; SZT2A
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001012961
Insert Size:	622 bp
Insert Sequence:	>SC207842 3'UTR clone of NM_001012961 The sequence shown below is from the reference sequence of NM_001012961. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATGTTGGCCAGGCTGCCCTCGAAGTTCTGACCTCAAGTGATCCACCTGCCTCAGCCTCCCAAAGTGCTG
GGATTACAGGCATGAGCTACTGTGCTGGCCAAAGCAGCATTGAACTGAGTTTGAAGAAAGAATTG
AGGTTAACCATTCAAAGAAGGGGAAATGGACTTTTCTGACAGAGGGACCAGCAAAAGCAAATACATGGT
ACTATGAACCTGTACTCTCAATCACAGGTGTTTATGGGAACCTACATATTATTCAGTGTCAACAGAGTAA
AATTTGAATGTGAAAGGGATGGTAAAGAGGAGACTGTAGACATAAGTGGGAGCTGGTTATGAAAGGCC
ATGAATGTGACACTAAGGACTTTGTGCTTTATCTATGGGGGAGTTCAAGTGTGGAGGATAAGAGCACTG
GATCTGGAATCACATTGCCTGTGTTTGCTTCTGCGCAGTACCCTTACTGGTTGTGTGCTTGTGACAA
GTTACTTAACCTTTCTGTGACCTGTGACCTTTTCCAAAGGGGATATTAACATTGCTTATATCATATG
GTGATTGAATTAATACATGTCAAGCACTTAGAAGTGTGCTTGCCAAAAAAAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001012961.1</u>
Locus ID:	149469
MW:	23.8