

Product datasheet for SC207853

CZIB (NM_017887) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CZIB (NM_017887) Human 3' UTR Clone
Symbol:	CZIB
Synonyms:	C1orf123
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017887
Insert Size:	590 bp
Insert Sequence:	>SC207853 3'UTR clone of NM_017887 The sequence shown below is from the reference sequence of NM_017887. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAGGTCACCCACCAGTTTGTGAAGTGCTGATCCCTCTTCCTTCCCAGTTGCCCTTAAGAACTGAGAAAG
GACAAAGTACTCTAAGCAGCAGAGCCACAGAGGCTCGTTCCTTTGACCCTTGCTCCTGGTGGCTATA
CGAAACCTTCACAATCTGCATGCTGGACTTTATTACAGCTTCCAAGCCCCATCAATAAAGCCCCCTGTT
CACGCTGCACTGGTGCATGAAGGTGAAGTAAGGAGGCAGCTACTCCCTCCACAAATAAGAATTCGTGCA
GCAGCAGTTCACCTTAGGAAAATGAAGAGTATTTTTAAGGAGATGGTATCTCCCCAGTCAAGGGCAA
GTTTACAGAGATGAAACCGCATCCACCATCACTCAATGAGCCAGAAGCCAATGGTGACGTACTGTTA
CCTGAACTTACTGAACTTGCTGTTGAAAGATCTTTGTGAGGATTAGGAATTAGGTTAAAAGAAATTAA
GAACCATCTTCAAGCAAAAATTAACCTTTATTTCTGCTTAAACAATAAATACACCTGAGTTAGTTTTCC
AAACCTTTCCTCCTGATTAATGCCCTTAAACTTAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_017887.3</u>
Locus ID:	54987
MW:	22.2