

## Product datasheet for **SC207856**

### SMN1 (NM\_022874) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	SMN1 (NM_022874) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SMN1
Synonyms:	BCD541; GEMIN1; SMA; SMA1; SMA2; SMA3; SMA4; SMA@; SMN; SMNT; T-BCD541; TDRD16A
ACCN:	NM_022874
Insert Size:	623 bp
Insert Sequence:	>SC207856 3'UTR clone of NM_022874

The sequence shown below is from the reference sequence of NM\_022874. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAGGAAGGTGCTCACATTCTTAAATTAAGGAGAAATGCTGGCATAGAGCAGCACTAAATGACACCAC
TAAAGAAACGATCAGACAGATCTGGAATGTGAAGCGTTATAGAAGATAACTGGCCTCATTCTTCAAAA
TATCAAGTGTGGGAAAGAAAAAGGAAGTGAATGGGTAACCTCTTCTTGATTAAGGTTATGTAATAA
CCAAATGCAATGTGAAATATTTTACTGGACTCTATTTTAAAAACCATCTGTAAGAACTGGGGTGGGG
GTGGGAGGCCAGCACGGTGGTGAGGCGAGTTGAGAAAAATTTGAATGTGGATTAGATTTTGAATGATATTG
GATAATTATTGGTAATTTTATGAGCTGTGAGAAGGGTGTGTAGTTTATAAAAGACTGTCTTAATTTGC
ATACTTAAGCATTAGGAATGAAGTGTAGAGTGTCTTAAATGTTTCAAATGGTTTAAACAAAATGTAT
GTGAGGCGTATGTGGCAAAATGTTACAGAATCTAACTGGTGGACATGGCTGTTTCATTGTACTGTTTTT
TCTATCTTCTATATGTTTAAAAAGTATATAAAAAATATTTAATTTTTTTTTTAAAAAAAAAAAAAAAA
AA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online >](#)

RefSeq: [NM\\_022874.2](#)

**Summary:** This gene is part of a 500 kb inverted duplication on chromosome 5q13. This duplicated region contains at least four genes and repetitive elements which make it prone to rearrangements and deletions. The repetitiveness and complexity of the sequence have also caused difficulty in determining the organization of this genomic region. The telomeric and centromeric copies of this gene are nearly identical and encode the same protein. However, mutations in this gene, the telomeric copy, are associated with spinal muscular atrophy; mutations in the centromeric copy do not lead to disease. The centromeric copy may be a modifier of disease caused by mutation in the telomeric copy. The critical sequence difference between the two genes is a single nucleotide in exon 7, which is thought to be an exon splice enhancer. Note that the nine exons of both the telomeric and centromeric copies are designated historically as exon 1, 2a, 2b, and 3-8. It is thought that gene conversion events may involve the two genes, leading to varying copy numbers of each gene. The protein encoded by this gene localizes to both the cytoplasm and the nucleus. Within the nucleus, the protein localizes to subnuclear bodies called gems which are found near coiled bodies containing high concentrations of small ribonucleoproteins (snRNPs). This protein forms heteromeric complexes with proteins such as SIP1 and GEMIN4, and also interacts with several proteins known to be involved in the biogenesis of snRNPs, such as hnRNP U protein and the small nucleolar RNA binding protein. Multiple transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2014]

Locus ID: 6606

MW: 24.1