

Product datasheet for **SC207906**

RAG1AP1 (SLC50A1) (NM_001122839) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RAG1AP1 (SLC50A1) (NM_001122839) Human 3' UTR Clone
Symbol:	RAG1AP1
Synonyms:	HsSWEET1; RAG1AP1; SCP; slv; SWEET1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001122839
Insert Size:	602 bp
Insert Sequence:	>SC207906 3'UTR clone of NM_001122839 The sequence shown below is from the reference sequence of NM_001122839. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GACAGGAACTACTGGCTCCTGCAAACCTAGGCTGCTCATCTGACCACTGGGCACCTTAGTGCCAACT
GAACCAAAGAGACCTCCTTGTTTCAGCTGGGCTGCTGTCCAGCTTCCAGGTGCAGTGGGTGTGGGA
ACAAGAGATGACTTTGAGGATAAAAGACCAAAGAAAAAGCTTTACTTAGATGATTGATTGGGCCTAG
GAGATGAAATCACTTTTTATTTTTTAGAGATTTTTTTTTTAATTTTGGAGGTTGGGTGCAATCTTTA
GAATATGCCTTAAAGGCCGGCGCGGTGGCTCACGCCTGTAATCCAGCACTTTGGGAGGCCAAGGTG
GGCGGATCGCCTGAGGTCAAGAGTTCAAGACCAACCTGACTAACATGGTGAACCCCATCTCTACTAAA
AATACAAAATTAGCCAGGCATGATGGCACATGCCTGTAATCCAGATACTTGGGAGGCTGAGGCAGGAG
AATTGCTTGAACCCAGGAGGTGGAGGTTGCAGTGAGCTGAGATCGTGCCATTGTGATATGAATATGCCT
TATATGCTGATATGAATATGCCTTAAATAAAGTGTTCCCAACCCCTGCC
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001122839.2
Summary:	Mediates sugar transport across membranes. May stimulate V(D)J recombination by the activation of RAG1.[UniProtKB/Swiss-Prot Function]
Locus ID:	55974
MW:	22.7