

Product datasheet for SC207923

C11orf17 (AKIP1) (NM_020642) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C11orf17 (AKIP1) (NM_020642) Human 3' UTR Clone
Symbol:	C11orf17
Synonyms:	BCA3; C11orf17
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020642
Insert Size:	629 bp
Insert Sequence:	>SC207923 3'UTR clone of NM_020642 The sequence shown below is from the reference sequence of NM_020642. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CAAAGCGTGGACCTGGTCTTCCCTGTGTGATGTTGACCATCACTGCCATCACATCACCTTTTTTAAAGT
AGTAAGAATAAAGCCACTGTATGATTCTCTTAATAGCTATACATTAATCCTGTTTTAGTGCTGACTGG
GTCAGCCTTCCGGGAAGTGGAGTCTGTCTTTTCAAGTGTCTTTTGTGTTGGTTGGTTTTTTTTTG
AGACAGTCTCACTCTGTTGCCAGGCTGGAGTGCAGTGGCGTGATCTCGGCTCACTGCAAGTTCCGCCT
CCCGGGTTCACACCATTCTCCTGCCTCAGCCTCCCGAGTAGCTGGCACTACAGGCACCCGCCACCATGC
CCGGCTATTTTTTTGTATTTTTAGTAGAGACGGGGTTTCAACATGTTGGCCAGGATGGTCTCGATCTC
TTGACCTCGTGATCCACCCACCTTGGCCTCCCAAAGTGTGGGATTACAGGCGTGAGCCACCGCGCCCG
GCCTCAGTGCCTTTTTAACTTGAGGGTGTAGAGGTCTCCACGCTTGTTCCTGAAAGTAATATAAT
GATGCTGTCTGAACAGGTTTACTGCTTGTCTTCCAAGTAAAGGTTAATTATGATAATAAAGAGATTGG
GAAATGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACGCGCCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_020642.4</u>
Summary:	This gene encodes a nuclear protein that interacts with protein kinase A catalytic subunit, and regulates the effect of the cAMP-dependent protein kinase signaling pathway on the NF-kappa-B activation cascade. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]
Locus ID:	56672
MW:	23.3