

Product datasheet for SC207988

ANKRD58 (SOWAHD) (NM_001105576) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ANKRD58 (SOWAHD) (NM_001105576) Human 3' UTR Clone
Symbol:	ANKRD58
Synonyms:	ANKRD58
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001105576
Insert Size:	617 bp
Insert Sequence:	>SC207988 3'UTR clone of NM_001105576 The sequence shown below is from the reference sequence of NM_001105576. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CATCTGTTCCCTCATTCCAGGACCGTTGACAGGGACAGAGACTGGAGAGCTAGGAGGGGCTGTGACAC
TGTGGCGATGGCTAGGTCCTGGGTTGTCCCGGTTCCACCGAAGGAGAGGCGCTTGGACGCTGCTTGG
GCCTGCAAGGAACAGAACACGTCGGGGTCCGACTCAGGTACTTGTCTCAGGTCTCCTGTAACCAACCGGC
CTGGAGGACCCGGGACTCGGGCACCACCTCACCAAGAGAGAGTGAAGGACCAAGCTGGCCTGGCTCCG
AGTTCCAAAGCTACAGGACTAAGGAGTTGGGAGCAGGGAGCGTGGTCTGCTTGGGAGAGGGCAAGTTA
AGCTTCCAGGGGCCATTTCTGGGCAGGCCGACGCGCTGGGTTTATTAGGAAACATTTCGCTAGAAGAATG
AGTTAAGATTGTAAACCACCAATGCAGAGAAAACGCCTAACTCTGCCGGCCTCGCTCGGCCATTAATGG
GTCTTGGGGTGCGGGTAGAGTCAGCCTCTGACAACCTCCTCTGAGACGACCCAGCCTTACTGGTACTT
TTTCTCATGTATCACAGTTACTTCTTATGTATATTAAGTGAATATGTGTTCTTTTCACATGA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001105576.3</u>
Locus ID:	347454
MW:	22.8