

Product datasheet for SC207997

LAYN (NM_178834) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: LAYN (NM_178834) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: LAYN
ACCN: NM_178834
Insert Size: 1317 bp
Insert Sequence: >SC207997 3'UTR clone of NM_178834

The sequence shown below is from the reference sequence of NM_178834. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
TGGGTGAAAAATGAAATATATGGTTATTAGGACATATAAAAACTGAAACTGACAACAATGAAAAAGAA
ATGATAAGCAAAATCCTCTTATTTTCTATAAGGAAAATACACAGAAGGTCTATGAACAAGCTTAGATCA
GGTCCTGTGGATGAGCATGTGGTCCCCACGACCTCCTGTTGGACCCCCAGTTTTGGCTGTATCCTTTA
TCCCAGCCAGTCATCCAGCTCGACCTTATGAGAAGGTACCTTGCCAGGTCTGGCACATAGTAGAGTCT
CAATAAATGTCACCTTGGTTGGTTGTATCTAACTTTTAAAGGGACAGAGCTTTACCTGGCAGTGATAAAGA
TGGGCTGTGGAGCTTGGAAAACCACTCTGTTTTCTTGTCTATACAGCAGCACATATTATCATACAG
ACAGAAAATCCAGAATCTTTTCAAAGCCCACATATGGTAGCACAGGCTGGCCTGTGCATCGGCAATTCT
CATATCTGTTTTTTCAAAGAATAATAAAATCAAATAAAGAGCAGGAAACAGAGTGTTAGTCTGTGTCT
ACAGCCCTTCTCTGCATGTGGCCACAGGGGACCTTTTTTTGTTTCTCCTGACATCCAGACTTGGAAAT
ATCTAACTACTTGCAAACTAAAAATGAGGCCAGGCCAGTGGCTGACGCCTGTAATCCCAGAACCTTG
GGAGACCAAGATTGGAGGATAGCTTGAGTTCAGGAGTTCAGACCTTCTGGGCAAAATAGTGAGACTC
TGACTCTACAAAAATTTAAAAATTAGCAGGGCATGGTGGCATGCGCCTGCAGTCCCAGCTACTCAGGA
GGCCGAGGTGGGGGAATTGCTTGAGCCCAGGAGGTTGAGGCTGCAGTGAGCTAGGATCACGCCACTGCT
CTCCTGCGTGGGCAACAGAATGAGACCTGTCTCTAAACAAAAATTTTTTTAAAGAAAATAAAAAATGA
AATTTGAGTCTCAGCTTCTAGGCTTATGTACTCCCATCCCTGATACTACACGTTGTTCTGTTTGCATGT
TTGTTACTGCAGCATCAACATAACATGACTAATAATGATGACAGTTCTACCTTTTGTGACAGACCTCA
TTGTTTCAATTGTGTTTGGCCTCAAATTCAGAGAAAAGTTTCTCAGCCCAGCATCAATCATATACACAT
TGGAGGTCTGAAAGCCATTATTTAAATCACCAGCAATTTCCCTCAATCAGCTTCATAATCTCATATTTT
AACCTTTCCTCATATGTCTATTTTTTAATAGCTCAGCAACTTTGCTTCACATATGGAGCCCCAACCTC
TAGAGT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_178834.5</u>
Summary:	Receptor for hyaluronate.[UniProtKB/Swiss-Prot Function]
Locus ID:	143903
MW:	50