

Product datasheet for **SC208070**

CD89 (FCAR) (NM_133272) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD89 (FCAR) (NM_133272) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	FCAR
Synonyms:	CD89; CTB-61M7.2; FcalphaRI
ACCN:	NM_133272
Insert Size:	1587 bp



[View online »](#)

Insert Sequence: >SC208070 3'UTR clone of NM_133272
The sequence shown below is from the reference sequence of NM_133272. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTTGCACGAACACCAAGTGTCTGCAAGTAACACCTGGAGGTGAAGGCAGAGAGGAGCCAGGACTGTGG
AGTCCGCAAAAGCTACTTGAAGGACACAAGAGAGAAAAGCTCACTAAGAAGCTTGAATCTACTTTTTTT
TTTTTTGAGACAGAGTCTGGCTCTGTCAACCAGGCTGGAGTGCAGTGGAGCAATCTCGGCTCATTGAA
CCTCTTGGGTTCAAGTATTCTTGTGCCTCAGCCTCCCAAGTAGCTGGAATTACAGGCACATACCACTG
CACCCAGCTAATTTTTGTATTTTAGTAGAGATGGGGTTTCACTGTGTTGGCCAGGCTGGTCTCGAACT
CCTGACCTCAGGTGATCCACCCACCTTGGCCTCCCAAAGTGTGAGATTATAGGCATGAGCCACCACGC
CTGGCCAGATGCATGTTCAAACCAATCAAATGGTGTCTTATGCAGGACTGATCGATTTGCCCCAC
CTTTCTGCACATAAGTTATGGTTTTCCATCTTATCTGTCTTCTGATTTTTATATCTGTTTAATTTCT
TCCTTCATTGTTCTTCTTTTTTTATTTATTTATTTATTTATTTATTTATTTATTTATTTGAGACAGAG
TCTCACTCTGTTGCCAGGCTGGAGTGCAGTGGCAGATCTCGGCTCACTGCAACCTCTGCCTCCTGGG
TTCAAGTGATTCTCCTGCCTCGGCCTCCCAAGTAGCTGGGATTGCAGGCTCCACCATCACGCCAGCT
ACTTTTACAGTATTTTGTAGTAGACGGGGTTTCACTCACATTGGCCAAGCTGGTCTCAAACCTTCTGACC
TCGTGATCTGCCCGCCTCGGCCTCCCAAAGTGTGGGATTACAGATGTGAGCCACTGCGCCAGCCTTC
TTTTTATATTTTAAATGTGTCTCCCAAATATAAATGGTTGGTAAGCATGCCAAATATATTCAATAA
CCCCCTCCTTTATTTTTTTTGTGAAGTGAGGCTCTCCCTATGTTGCCTAAGCTGGTCTTGAAGTCC
TGGTCTCAAGCAATCCTCCTACCTCAGCCTCCTGTGTTCATCTACAAATTGATAAGAGTGAAGTC
ATAATCCTACAGGAGGATTACCCTATTTATTTCAAAACCTATTTCTACCGGATTTTCATACAAGGAA
TACAGGCATGTGTTTACCTCATTAAATTTATTTTTCACTTAGTTTTGATGATATTCACATATATTATC
AAGTGTGCAAAACATTAATTTCTTGTGTACAAAACCTCAAATGGTCTTCAAATAAATCCCCATTCTTTTT
TCTTATAAACTTTCACAGCTTTACCCTTGACAGACTTTACTCAAGGAAATCTAAGTTGGTCATATGTGG
CTCTTTCAGTATTGCTATTTACTTTCATTGTCCAGTAGCTTATGTATGAAAATATAAATTAAAAATGTA
AGGGTCTACTTCCAGTGAAGTGAAGGGACTTAGGCCACTTTTATCCTTTACTGAGAGCTTATCTCT
ACTTGATAAAATTTCTACTGTATTCTTGGCTTAACTCAGGTCCTGTGATTAATAAAAAAAAAATGCAAAGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_133272.4](#)

Summary:

This gene is a member of the immunoglobulin gene superfamily and encodes a receptor for the Fc region of IgA. The receptor is a transmembrane glycoprotein present on the surface of myeloid lineage cells such as neutrophils, monocytes, macrophages, and eosinophils, where it mediates immunologic responses to pathogens. It interacts with IgA-opsionized targets and triggers several immunologic defense processes, including phagocytosis, antibody-dependent cell-mediated cytotoxicity, and stimulation of the release of inflammatory mediators. Multiple alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Locus ID:

2204

MW:

61.5