

Product datasheet for **SC208073**

CD89 (FCAR) (NM_133269) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CD89 (FCAR) (NM_133269) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	FCAR
Synonyms:	CD89; CTB-61M7.2; FcalphaRI
ACCN:	NM_133269
Insert Size:	1587 bp



[View online »](#)

Insert Sequence: >SC208073 3'UTR clone of NM_133269
 The sequence shown below is from the reference sequence of NM_133269. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTTGCACGAACACCAAGTGTCTGCAAGTAACACCTGGAGGTGAAGGCAGAGAGGAGCCAGGACTGTGG
AGTCCGCAAAAGCTACTTGAAGGACACAAGAGAGAAAAGCTCACTAAGAAGCTTGAATCTACTTTTTTT
TTTTTTGAGACAGAGTCTGGCTCTGTCAACCAGGCTGGAGTGCAGTGGAGCAATCTCGGCTCATTGAA
CCTCTTGGGTTCAAGTATTCTTGTGCCTCAGCCTCCCAAGTAGCTGGAATTACAGGCACATACCACTG
CACCCAGCTAATTTTTGTATTTTTAGTAGAGATGGGGTTTCACTGTGTTGGCCAGGCTGGTCTCGAACT
CCTGACCTCAGGTGATCCACCCACCTTGGCCTCCCAAAGTGTGAGATTATAGGCATGAGCCACCACGC
CTGGCCAGATGCATGTTCAAACCAATCAAATGGTGTCTTATGCAGGACTGATCGATTTGCCCCAC
CTTTCTGCACATAAGTTATGGTTTTCCATCTTATCTGTCTTCTGATTTTTATATCTGTTTAATTTCT
TCCTTCATTGTTCTTCTTTTTTTATTTATTTATTTATTTATTTATTTATTTATTTATTTGAGACAGAG
TCTCACTCTGTTGCCAGGCTGGAGTGCAGTGGCAGATCTCGGCTCACTGCAACCTCTGCCTCCTGGG
TTCAAGTGATTCTCCTGCCTCGGCCTCCCAAGTAGCTGGGATTGCAGGCTCCACCATCACGCCAGCT
ACTTTTACAGTATTTTTAGTAGAGACGGGGTTTCATCACATTGGCCAAGCTGGTCTCAAACCTTCTGACC
TCGTGATCTGCCCGCCTCGGCCTCCCAAAGTGTGGGATTACAGATGTGAGCCACTGCGCCAGCCTTC
TTTTTATATTTTTAAATGTGTCTCCCAAATATAAATGGTTGGTAAGCATGCCAAATATATTCAATAA
CCCCCTCCTTTATTTTTTTTTGTTGAAGTGAGGCTCTCCCTATGTTGCCTAAGCTGGTCTTGAACCTC
TGGTCTCAAGCAATCCTCCTACCTCAGCCTCCTGTGTTCATCTACAAATTGATAAGAGTGAAGTC
ATAATCCTACAGGAGGATTACCCTATTTATTTACAAACCTATTTCTACCGGATTTTCATACAAGGAA
TACAGGCATGTGTTTACCTCATTAAATTTATTTTTCACTTAGTTTTGATGATATTCACATATATTATC
AAGTGTGCAAAACATTAATTTCTTGTGTACAAAACCTCAAATGGTCTTCAAATAAATCCCCATTCTTTTT
TCTTATAAACTTTACAGCTTTACCCTTGACAGACTTTACTCAAGGAAATCTAAGTTGGTCATATGTGG
CTCTTTCAGTATTGCTATTTACTTTCATTGTCCAGTAGCTTATGTATGAAAATATAAATTAAAAATGTA
AGGGTCTACTTCCAGTGAAGTGAAGGGACTTAGGCCACTTTTATCCTTTACTGAGAGCTTATCTCT
ACTTGATAAAATTTCTACTGTATTCTTGGCTTAACTCAGGTCCTGTGATTAATAAAAAAAAAATGCAAAGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_133269.4](#)

Summary:

This gene is a member of the immunoglobulin gene superfamily and encodes a receptor for the Fc region of IgA. The receptor is a transmembrane glycoprotein present on the surface of myeloid lineage cells such as neutrophils, monocytes, macrophages, and eosinophils, where it mediates immunologic responses to pathogens. It interacts with IgA-opsonized targets and triggers several immunologic defense processes, including phagocytosis, antibody-dependent cell-mediated cytotoxicity, and stimulation of the release of inflammatory mediators. Multiple alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Locus ID:

2204

MW:

61.5