

Product datasheet for SC208122

FILIP1 (NM_015687) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FILIP1 (NM_015687) Human 3' UTR Clone
Symbol:	FILIP1
Synonyms:	FILIP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015687
Insert Size:	631 bp
Insert Sequence:	>SC208122 3'UTR clone of NM_015687 The sequence shown below is from the reference sequence of NM_015687. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ACCACTCTCTCGGAGGGGGGAAGGGCTAGGGCAGTGGCTAAGGGGGTATGTTGTGCAGATGCTACTG
CTGCCGTGAAAGTGAACCTTCATCTGTTTGTGCCAGTCTTTACATGTAATAATTAAGTTTAAATAT
TGTGTTTATAAAATAACCAACTAATAACCATTGTCTTTCCATTTGTGCATTTGTTTGTATGCTGGG
GAACAAAATTAGCAAACTATTGCTTGCTGCCTAGAAGCCAGGGCGTGGTTTCTAGTTCCAGTTTGTCT
TCTAGCAAGTGGACCCATCAATAGACCCATCTGAGCCTGTTTCTCATCAGTTAGATGTGGGGACTCAA
TCACACGCTCTTCAAGTCCGGCTCCCATATTTCTAATTGCAAGCCAAATTTAATGTACCTTGTTCAC
AATAATTTTTTATTAATAAATCCTATTACAAAATAAGACATACTTTAACTATTGTCAATTGCTCTTT
CACATCATGAATTTGCTTTATGTGCTGGAAAAACATCAGCTATCAGGGCCTGGACCTCTAAA
ATTTTGCAAAAACAAAAGTTCTAAGATGATTTAGGAAATAATGTGAACATGTAATAAATGGAATG
AAATATGGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_015687.5
Summary:	This gene encodes a filamin A binding protein. The encoded protein promotes the degradation of filamin A and may regulate cortical neuron migration and dendritic spine morphology. Mice lacking a functional copy of this gene exhibit reduced dendritic spine length and altered excitatory signaling. [provided by RefSeq, Oct 2016]
Locus ID:	27145
MW:	24.4