

Product datasheet for **SC208128**

KDM3A / JHDM2A (KDM3A) (NM_018433) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KDM3A / JHDM2A (KDM3A) (NM_018433) Human 3' UTR Clone
Symbol:	KDM3A / JHDM2A
Synonyms:	JHDM2A; JHMD2A; JMJD1; JMJD1A; TSGA
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018433
Insert Size:	627 bp
Insert Sequence:	>SC208128 3'UTR clone of NM_018433 The sequence shown below is from the reference sequence of NM_018433. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GCCAGTGAATCCAGTTTTGGCAAACCTTAATCTCCCTGCACATTGGAAATGAATTACAGGCAGCTGTTC
AAACTCTTCAGGCAGGATTCTGTGGACTTTGAGATTCATGTTACCTCATCTTCTTTTAACTGTAC
CCAACCTGTGAGGGTACTCTGTCTAATGTATATTTCTAGTGTTCACAGACAGTAAATGTGTATATGTAG
TAACTATTTACAGAACATGCATCCTTAACTGTGACTTCTCACCTAGTGCAGAACTTTACCAGGCTGT
AAAAGCAAAACCTCGTATCAGCTCTGGAACAATACCTGCAGTTATTCTTCAGCTGTTTGGACAACCTAG
ATTGGGTTTATACTATTAGGAATCACTGCACAGTTTATTTGGGTTGTGTTTGTGTCTGAGTCCCTC
CCTCATCCCTTAGGGTCCAGAAGAGCAATGGAGGAAGTGACAGCTAATGTTGCAGTTCTTATTGTATGG
CATAGGACTGGCATTATATAGCAGAAATCAACTACTGTACAATTTCTTGGGGTTAACCATCTTTAGTTA
AATGGAATTTAATTTAAATGACGCTTTGCTAATTTAAGTGTTAAGCATTTTGCATTAAATATTCAAT
ATAATA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018433.6
Summary:	This gene encodes a zinc finger protein that contains a jumonji domain and may play a role in hormone-dependent transcriptional activation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]
Locus ID:	55818
MW:	23.6