

Product datasheet for **SC208140**

DIS3L (NM_133375) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	DIS3L
Synonyms:	DIS3L1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_133375
Insert Size:	616 bp
Insert Sequence:	>SC208140 3'UTR clone of NM_133375 The sequence shown below is from the reference sequence of NM_133375. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CTGGATGTTTCAAACAATTATGGAATAAGAGGCTCTTACTTCACTAAGAGCTGTCATATGTGAATGT TTTACAGTCTTTTCAAACCTAACATTTAATGTGTGTCAGTCTAGTCGATCAGGACTGGGTAG CTATTTTCGCATATATGTAATGTTCTCAGCCGGGCACGGTGGCTCAGCCTGTAACCCAGCACTTTG GGAGGCTGAGGCGGGCGGATCACGAGGTCAGGAGATTGAGACCATCTGGCTAACACGGTGAAACCCAG TCTCTACTAAAAATACAAAAATTAGCCCGCGTGGTGGCATGCGCCTGTAGTCCAGCTACTTGGGAGG CTGAAGCAGGAGAATTGCCTGAACCCAGGAAGAGGAGTTGCAGTGAGCCGAGATCGCACCCTGCATT CCAGCCTGGGCAACAGAGCGAGACTCCATCTCAAAAAAAAAAAAAAAAAAAGTTCTCTCATTTC ATTAAAGTTGCATTAAATAAAGTATAATTAGGTCACTATGGAACAGAGTTTTCAGTAATGAGTGGACA GTAAGTGGTGGCTCTGCAGATGGCCCTTTCTAATAAGTTTAATAAACCCGAAATTACTGGTTA ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_133375.5</u>
Summary:	The cytoplasmic RNA exosome complex degrades unstable mRNAs and is involved in the regular turnover of other mRNAs. The protein encoded by this gene contains 3'-5' exoribonuclease activity and is a catalytic component of this complex. [provided by RefSeq, May 2016]
Locus ID:	115752
MW:	22.4